

BRD3 shRNA Plasmid (h): sc-60284-SH

BACKGROUND

The bromodomain-containing proteins include BRD2, BRD3, BRD4 and BRDT. BRD2 (RING3 protein) is a mitogen-activated nuclear protein whose gene is located in the human MHC II region, suggesting its relation to HLA-associated diseases. The gene encoding BRD3 (RING3-like protein) contains two bromodomains and maps to chromosome 9q34.2. BRD4 (HUNK1 protein) is a nuclear protein involved in the regulation of chromosomal dynamics during mitosis. The testis-specific bromodomain protein BRDT contains a PEST sequence, indicating that it undergoes rapid intracellular degradation. The bromodomain-containing proteins are ubiquitously expressed.

REFERENCES

1. Thorpe, K.L., et al. 1997. Chromosomal localization, gene structure and transcription pattern of the ORFX gene, a homologue of the MHC-linked RING3 gene. *Gene* 200: 177-183.
2. Zhou, M., et al. 2003. Expression of BRD7-interacting proteins, BRD2 and BRD3, in nasopharyngeal carcinoma tissues. *Ai Zheng* 22: 123-127.
3. Shang, E., et al. 2004. Identification of unique, differentiation stage-specific patterns of expression of the bromodomain-containing genes Brd2, Brd3, Brd4, and Brdt in the mouse testis. *Gene Expr. Patterns* 4: 513-519.
4. Boyer, A., et al. 2004. Pre-sertoli specific gene expression profiling reveals differential expression of Ppt1 and Brd3 genes within the mouse genital ridge at the time of sex determination. *Biol. Reprod.* 71: 820-827.
5. Crowley, T., et al. 2004. Change in nuclear-cytoplasmic localization of a double-bromodomain protein during proliferation and differentiation of mouse spinal cord and dorsal root ganglia. *Brain Res. Dev. Brain Res.* 149: 93-101.

CHROMOSOMAL LOCATION

Genetic locus: BRD3 (human) mapping to 9q34.2.

PRODUCT

BRD3 shRNA Plasmid (h) is a pool of 3 target-specific lentiviral vector plasmids each encoding 19-25 nt (plus hairpin) shRNAs designed to knock down gene expression. Each plasmid contains a puromycin resistance gene for the selection of cells stably expressing shRNA. Each vial contains 20 µg of lyophilized shRNA plasmid DNA. Suitable for up to 20 transfections. Also see BRD3 siRNA (h): sc-60284 and BRD3 shRNA (h) Lentiviral Particles: sc-60284-V as alternate gene silencing products.

RESEARCH USE

The purchase of this product conveys to the buyer the nontransferable right to use the purchased amount of the product and all replicates and derivatives for research purposes conducted by the buyer in his laboratory only (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party, or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes.

STORAGE AND RESUSPENSION

Store lyophilized shRNA plasmid DNA at 4° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at 4° C for short term storage or -80° C for long term storage. Avoid repeated freeze thaw cycles.

Resuspend lyophilized shRNA plasmid DNA in 200 µl of the deionized water provided. Resuspension of the shRNA plasmid DNA in 200 µl of deionized water makes a 0.1 µg/µl solution in a 10 mM Tris, 1 mM EDTA buffered solution.

APPLICATIONS

BRD3 shRNA Plasmid (h) is recommended for the inhibition of BRD3 expression in human cells.

SUPPORT REAGENTS

For optimal shRNA Plasmid transfection efficiency, Santa Cruz Biotechnology's shRNA Plasmid Transfection Reagent: sc-108061 (0.2 ml) and shRNA Plasmid Transfection Medium: sc-108062 (20 ml) are recommended. Control shRNAs are available as 20 µg lyophilized plasmid DNA. Each encodes a scrambled shRNA sequence that will not lead to the specific degradation of any known cellular mRNA. Control shRNA Plasmids include: sc-108060, sc-108065 and sc-108066.

GENE EXPRESSION MONITORING

BRD3 (B-12): sc-515666 is recommended as a control antibody for monitoring of BRD3 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgGκ BP-FITC: sc-516140 or m-IgGκ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor BRD3 gene expression knockdown using RT-PCR Primer: BRD3 (h)-PR: sc-60284-PR (20 µl). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.