

CARD 6 siRNA (h): sc-60332

BACKGROUND

Membrane-associated guanylate kinase (MAGUK) family members localize to the plasma membrane and function as molecular scaffolds for the assembly of multi-protein complexes. The MAGUK family includes several mammalian proteins related to the *Drosophila* tumor suppressor discs-large (dlg) gene product, such as postsynaptic proteins, GKAPs, the tight junction associated proteins (ZO-1-3) and the caspase-associated recruitment domain (CARD) proteins (CARD 6, CARD 8-12 and CARD 14). CARD 6 is a 311 amino acid protein that associates with microtubules and acts as a selective mediator of NF κ B activation. CARD 6 is structurally and functionally related to the super-family of interferon (IFN)-inducible GTPases, a group that comprises host defense proteins which mediate cell-autonomous immunity. RICK, the receptor-interacting protein (RIP)-like interacting caspase-like apoptosis regulatory protein kinase, interacts with CARD 6 and targets it to aggresomes, regions of the cell where protein aggregates collect.

REFERENCES

1. Online Mendelian Inheritance in Man, OMIM[™]. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 609986. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
2. Stehlik, C., et al. 2003. CARD 6 is a modulator of NF κ B activation pathways. *J. Biol. Chem.* 278: 31941-31949.
3. Carson, S.D. 2004. Coxsackievirus and adenovirus receptor (CAR) is modified and shed in membrane vesicles. *Biochemistry* 43: 8136-8142.
4. Hong, G.S., et al. 2005. Caspase recruitment domain (CARD) as a bi-functional switch of caspase regulation and NF κ B signals. *J. Biochem. Mol. Biol.* 35: 19-23.
5. Dufner, A., et al. 2006. Caspase recruitment domain protein 6 is a microtubule-interacting protein that positively modulates NF κ B activation. *Proc. Natl. Acad. Sci. USA* 103: 988-993.
6. Chen, Y.R., et al. 2006. Substitutions of prolines examine their role in kinetic trap formation of the caspase recruitment domain (CARD) of RICK. *Protein Sci.* 15: 395-409.

CHROMOSOMAL LOCATION

Genetic locus: CARD6 (human) mapping to 5p13.1.

PRODUCT

CARD 6 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see CARD 6 shRNA Plasmid (h): sc-60332-SH and CARD 6 shRNA (h) Lentiviral Particles: sc-60332-V as alternate gene silencing products.

For independent verification of CARD 6 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-60332A, sc-60332B and sc-60332C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

CARD 6 siRNA (h) is recommended for the inhibition of CARD 6 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor CARD 6 gene expression knockdown using RT-PCR Primer: CARD 6 (h)-PR: sc-60332-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.