



# CLEC-6 siRNA (h): sc-60398

## BACKGROUND

The human  $\beta$ -glucan protein (Dectin-1) is a small, type II transmembrane receptor that enables  $\beta$ -glucan dependent, nonopsonic recognition of zymosan and other yeast-derived particles by primary macrophages. Dectin-1 is expressed in dendritic cells and is the human homolog of the C-type (calcium dependent) lectin-like receptor (CLEC) family that plays an important role in regulating innate immunity. The CLEC protein structure has a specific fold that provides a highly static scaffold for combinatorial display of variable antigen residues. This fold differs from the classic immunoglobulin fold, illustrating an evolutionary solution for balancing diversity against stability. CLEC-6 is a single-pass, type II transmembrane protein that is highly expressed in bone marrow and spleen and weakly expressed in peripheral blood leukocytes. CLEC-6 may function as a receptor in the antigen uptake at the infection location.

## REFERENCES

1. Arce, I., et al. 2001. Molecular and genomic characterization of human DLEC, a novel member of the C-type lectin receptor gene family preferentially expressed on monocyte-derived dendritic cells. *Eur. J. Immunol.* 31: 2733-2740.
2. Hermanz-Falcón, P., et al. 2001. Cloning of human Dectin-1, a novel C-type lectin-like receptor gene expressed on dendritic cells. *Immunogenetics* 53: 288-295.
3. Sobanov, Y., et al. 2001. A novel cluster of lectin-like receptor genes expressed in monocytic, dendritic and endothelial cells maps close to the NK receptor genes in the human NK gene complex. *Eur. J. Immunol.* 31: 3493-3503.
4. Yokota, K., et al. 2001. Identification of a human homologue of the dendritic cell-associated C-type lectin-1, Dectin-1. *Gene* 272: 51-60.
5. Gavino, A.C., et al. 2005. Identification and expression profiling of a human C-type lectin, structurally homologous to mouse Dectin-2. *Exp. Dermatol.* 14: 281-288.
6. McMahon, S.A., et al. 2005. The C-type lectin fold as an evolutionary solution for massive sequence variation. *Nat. Struct. Mol. Biol.* 12: 886-892.
7. Chaipan, C., et al. 2006. DC-SIGN and CLEC-2 mediate human immunodeficiency virus type 1 capture by platelets. *J. Virol.* 80: 8951-8960.
8. Watson, A.A. and O'Callaghan, C.A. 2006. Crystallization and X-ray diffraction analysis of human CLEC-2. *Acta Crystallogr. Sect. F, Struct. Biol. Cryst. Commun.* 61: 1094-1096.

## CHROMOSOMAL LOCATION

Genetic locus: CLEC4D (human) mapping to 12p13.31.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

## PRODUCT

CLEC-6 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see CLEC-6 shRNA Plasmid (h): sc-60398-SH and CLEC-6 shRNA (h) Lentiviral Particles: sc-60398-V as alternate gene silencing products.

For independent verification of CLEC-6 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-60398A, sc-60398B and sc-60398C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

CLEC-6 siRNA (h) is recommended for the inhibition of CLEC-6 expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor CLEC-6 gene expression knockdown using RT-PCR Primer: CLEC-6 (h)-PR: sc-60398-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.