

Olfr558 siRNA (m): sc-60504

BACKGROUND

G protein-coupled receptors (GPCRs), also designated seven transmembrane (7TM) receptors and heptahelical receptors, are a protein family which interact with G proteins (heterotrimeric GTPases) to synthesize intracellular second messengers such as diacylglycerol, cyclic AMP, inositol phosphates, and calcium ions. Their diverse biological functions range from vision and olfaction to neuronal and endocrine signaling and are involved in many pathological conditions. The GPCR family represents the largest class of targets for modern drugs. Dresden G protein-coupled receptor (D-GPCR) belongs to the subfamily of odorant-like orphan G protein-coupled receptors. D-GRCP is selectively overexpressed in human prostate cancer (PCA). The expression of D-GRCP increases with higher tumor grades and stages. D-GRCP is a potential tumor marker and a target for prostate cancer therapy.

REFERENCES

1. Lamah, J., Cone, R.I., Maeda, S., Philip, M., Corbani, M., Nádasdi, L., Ramachandran, J., Smith, G.M. and Sadée, W. 1991. Structure and function of G protein-coupled receptors. *Pharm. Res.* 7: 1213-1221.
2. Probst, W.C., Snyder, L.A., Schuster, D.I., Brosius, J. and Sealfon, S.C. 1992. Sequence alignment of the G protein-coupled receptor superfamily. *DNA Cell Biol.* 11: 1-20.
3. George, S.R., O'Dowd, B.F. and Lee, S.P. 2002. G protein-coupled receptor oligomerization and its potential for drug discovery. *Nat. Rev. Drug Discov.* 1: 808-820.
4. Harmor, A.J. 2002. Family-B G protein-coupled receptors. *Genome Biol.* 2: REVIEWS3013.
5. Weigle, B., Fuessel, S., Ebner, R., Temme, A., Schmitz, M., Schwind, S., Kiessling, A., Rieger, M.A., Meye, A., Bachmann, M., Wirth, M.P. and Rieber, E.P. 2004. D-GPCR: a novel putative G protein-coupled receptor overexpressed in prostate cancer and prostate. *Biochem. Biophys. Res. Commun.* 322: 239-249.
6. Fuessel, S., Weigle, B., Schmidt, U., Baretton, G., Koch, R., Bachmann, M., Rieber, E.P., Wirth, M.P. and Meye, A. 2006. Transcript quantification of Dresden G protein-coupled receptor (D-GPCR) in primary prostate cancer tissue pairs. *Cancer Lett.* 236: 95-104.

CHROMOSOMAL LOCATION

Genetic locus: Olfr558 (mouse) mapping to 7 E3.

PRODUCT

Olfr558 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Olfr558 shRNA Plasmid (m): sc-60504-SH and Olfr558 shRNA (m) Lentiviral Particles: sc-60504-V as alternate gene silencing products.

For independent verification of Olfr558 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-60504A, sc-60504B and sc-60504C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Olfr558 siRNA (m) is recommended for the inhibition of Olfr558 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Olfr558 gene expression knockdown using RT-PCR Primer: Olfr558 (m)-PR: sc-60504-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.