

Pancreatic Lipase siRNA (h): sc-61285

BACKGROUND

The lipase gene family belongs to one of the most robust genetic superfamilies found in living organisms, which includes esterases and thioesterases. Members of the AB hydrolase subfamily include hepatic lipase (HL), endothelial lipase (EL), lipoprotein lipase (LPL), Pancreatic Lipase (PL), gastric lipase (GL) and LCAT. These family members play a crucial role in the metabolism of lipids. Pancreatic lipase, also designated pancreatic triacylglycerol acyl hydrolase, is important for dietary fat absorption as it hydrolyses triglycerides into diglycerides, monoglycerides and free fatty acids.

REFERENCES

1. Lowe, M.E., et al. 1989. Cloning and characterization of human Pancreatic Lipase cDNA. *J. Biol. Chem.* 264: 20042-20048.
2. Winkler, F.K., et al. 1990. Structure of human Pancreatic Lipase. *Nature* 343: 771-774.
3. Yajima, H., et al. 2005. Prevention of diet-induced obesity by dietary isomerized hop extract containing isohumulones, in rodents. *Int. J. Obes.* 29: 991-997.
4. Han, L.K., et al. 2005. Anti-obesity effects of chikusetsusaponins isolated from *Panax japonicus* rhizomes. *BMC Complement. Altern. Med.* 5: 9.
5. Bijvelds, M.J., et al. 2005. Fat absorption in cystic fibrosis mice is impeded by defective lipolysis and post-lipolytic events. *Am. J. Physiol. Gastrointest. Liver Physiol.* 288: G646-G653.
6. Sharma, N., et al. 2005. Screening of some medicinal plants for anti-lipase activity. *J. Ethnopharmacol.* 97: 453-456.
7. Pappan, K.L., et al. 2005. Pancreatic β -cell lipoprotein lipase independently regulates islet glucose metabolism and normal Insulin secretion. *J. Biol. Chem.* 280: 9023-9029.

CHROMOSOMAL LOCATION

Genetic locus: PNLIP (human) mapping to 10q25.3.

PRODUCT

Pancreatic Lipase siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Pancreatic Lipase shRNA Plasmid (h): sc-61285-SH and Pancreatic Lipase shRNA (h) Lentiviral Particles: sc-61285-V as alternate gene silencing products.

For independent verification of Pancreatic Lipase (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-61285A, sc-61285B and sc-61285C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Pancreatic Lipase siRNA (h) is recommended for the inhibition of Pancreatic Lipase expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

Pancreatic Lipase (A-3): sc-374612 is recommended as a control antibody for monitoring of Pancreatic Lipase gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Pancreatic Lipase gene expression knockdown using RT-PCR Primer: Pancreatic Lipase (h)-PR: sc-61285-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.