

# Pepsin C siRNA (m): sc-61319

## BACKGROUND

Pepsin is one of the main proteolytic enzymes secreted by the gastric mucosa. Pepsin consists of a single polypeptide chain and arises from its precursor, pepsinogen, by removal of a 41 amino acid segment from the N-terminus. Pepsinogen is synthesized in the stomach lining, and hydrochloric acid, also produced by the gastric mucosa, is necessary to convert the inactive enzyme and to maintain the optimum acidity (pH 1-3) for pepsin function. Pepsin is particularly effective in cleaving peptide bonds involving aromatic amino acids. Pepsin shows extremely broad specificity, and although bonds involving phenylalanine and leucine are preferred, many others are also cleaved to some extent. The amino acid composition of Pepsin C differs from those of pepsinogen and pepsin especially in the content of basic amino acids, glutamic acid, aspartic acid, leucine and isoleucine.

## REFERENCES

1. Ryle, A.P. and Hamilton, M.P. 1967. Pepsinogen C and Pepsin C. Further purification and amino acid composition. *Biochem. J.* 101: 176-183.
2. Ryle, A.P., Leclerc, J. and Falla, F. 1969. The substrate specificity of pepsin C. *Biochem. J.* 110: 4P.
3. Kay, J. and Ryle, A.P. 1972. An active site peptide from pepsin C. *Biochem. J.* 123: 75-82.
4. Kageyama, T. and Takahashi, K. 1977. Pepsinogen C and pepsin C from gastric mucosa of Japanese monkey. Purification and characterization. *J. Biochem.* 80: 983-992.
5. Auffret, C.A. and Ryle, A.P. 1979. The catalytic activity of pig pepsin C towards small synthetic substrates. *Biochem. J.* 179: 239-246.
6. Szecsi, P.B., Dalgaard, D., Stakemann, G., Wagner, G. and Foltmann, B. 1989. The concentration of pepsinogen C in human semen and the physiological activation of zymogen in the vagina. *Biol. Reprod.* 40: 653-659.
7. Sánchez, L.M., Freije, J.P., Merino, A.M., Vizoso, F., Foltmann, B. and López-Otín, C. 1992. Isolation and characterization of a pepsin C zymogen produced by human breast tissues. *J. Biol. Chem.* 267: 24725-24731.
8. Schreiber, S., Buckner, R., Groll, C., Azevedo-Vethacke, M., Scheid, P., Gattermann, S., Josenhans, C. and Suerbaum, S. 2006. Gastric antibacterial efficiency is different for pepsin A and C. *Arch. Microbiol.* 184: 335-340.

## CHROMOSOMAL LOCATION

Genetic locus: Pgc (mouse) mapping to 17 C.

## PRODUCT

Pepsin C siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Pepsin C shRNA Plasmid (m): sc-61319-SH and Pepsin C shRNA (m) Lentiviral Particles: sc-61319-V as alternate gene silencing products.

For independent verification of Pepsin C (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-61319A, sc-61319B and sc-61319C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

Pepsin C siRNA (m) is recommended for the inhibition of Pepsin C expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

Pepsin C (E-9): sc-374044 is recommended as a control antibody for monitoring of Pepsin C gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Pepsin C gene expression knockdown using RT-PCR Primer: Pepsin C (m)-PR: sc-61319-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.