



RGS20 siRNA (h): sc-61470

BACKGROUND

The regulators of G protein signaling (RGS) proteins inhibit heterotrimeric G protein signaling. RGS proteins work by functioning as GTPase-activating proteins (which increase the GTPase activity of G protein α subunits) thereby driving G proteins into their inactive GDP-bound form. RGS20 is expressed exclusively in brain, with highest levels in the caudate nucleus and temporal lobe. RGS20 belongs to the RZ subfamily because it is highly selective for the α subunit on G proteins. However, if protein kinase C phosphorylates the α subunit, the G protein is much less susceptible to RGS20 action. RGS20 directly interacts with the microtubule-destabilizing protein SCG10 (superior cervical ganglia, neural specific 10) and blocks its ability to induce microtubule disassembly.

REFERENCES

1. Glick, J.L., Meigs, T.E., Miron, A. and Casey, P.J. 1998. RGSZ1, a G_z -selective regulator of G protein signaling whose action is sensitive to the phosphorylation state of G_{α_z} . *J. Biol. Chem.* 273: 26008-26013.
2. Wang, J., Ducret, A., Tu, Y., Kozasa, T., Aebersold, R. and Ross, E.M. 1998. RGSZ1, a G_z -selective RGS protein in brain. Structure, membrane association, regulation by G_{α_z} phosphorylation, and relationship to a G_z gtpase-activating protein subfamily. *J. Biol. Chem.* 273: 26014-26025.
3. Barker, S.A., Wang, J., Sierra, D.A. and Ross, E.M. 2001. RGSZ1 and Ret RGS: two of several splice variants from the gene RGS20. *Genomics* 78: 223-229.
4. Nixon, A.B., Grenningloh, G. and Casey, P.J. 2002. The interaction of RGSZ1 with SCG10 attenuates the ability of SCG10 to promote microtubule disassembly. *J. Biol. Chem.* 277: 18127-18133.
5. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 607193. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Wang, Y. and Young, K.H. 2004. Analysis of RGSZ1 protein interaction with G_{α_i} subunits. *Methods Enzymol.* 390: 31-52.
7. Nixon, A.B. and Casey, P.J. 2004. Analysis of the regulation of microtubule dynamics by interaction of RGSZ1 (RGS20) with the neuronal stathmin, SCG10. *Methods Enzymol.* 390: 53-64.
8. Garzón, J., Rodríguez-Muñoz, M., López-Fando, A. and Sánchez-Blázquez, P. 2005. The RGSZ2 protein exists in a complex with μ -opioid receptors and regulates the desensitizing capacity of G_z proteins. *Neuropsychopharmacology* 30: 1632-1648.

CHROMOSOMAL LOCATION

Genetic locus: RGS20 (human) mapping to 8q11.23.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

RGS20 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see RGS20 shRNA Plasmid (h): sc-61470-SH and RGS20 shRNA (h) Lentiviral Particles: sc-61470-V as alternate gene silencing products.

For independent verification of RGS20 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-61470A, sc-61470B and sc-61470C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

RGS20 siRNA (h) is recommended for the inhibition of RGS20 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor RGS20 gene expression knockdown using RT-PCR Primer: RGS20 (h)-PR: sc-61470-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.