

SCO1 siRNA (h): sc-61505

BACKGROUND

The SCO1 and SCO2 protein homologs belong to the SCO1/2 family of proteins. SCO1 and SCO2 both localize to the mitochondrion and are inner membrane proteins crucial for copper insertion or transport to the active site of cytochrome c oxidase (COX). COX is a crucial component in energy production because it functions as the terminal enzyme in the respiratory chain. SCO1 is predominantly expressed in highly oxidative phosphorylation tissues such as brain, heart and muscle, while SCO2 is ubiquitously expressed. Defects in the gene encoding for SCO1 may cause cytochrome c oxidase deficiency, a hetero-genous disorder. Defects in the gene encoding for SCO2 may cause cardioencephalomyopathy with cytochrome c oxidase deficiency, a fatal infantile disorder characterized by hypertrophic cardiomyopathy, lactic acidosis and gliosis.

REFERENCES

1. Jaksch, M., et al. 2001. Cytochrome c oxidase deficiency due to mutations in SCO2, encoding a mitochondrial copper-binding protein, is rescued by copper in human myoblasts. *Hum. Mol. Genet.* 10: 3025-3035.
2. Balatri, E., et al. 2003. Solution structure of SCO1: a thioredoxin-like protein involved in cytochrome c oxidase assembly. *Structure* 11: 1431-1443.
3. Horng, Y.C., et al. 2004. Specific copper transfer from the COX17 metallo-chaperone to both SCO1 and COX11 in the assembly of yeast cytochrome c oxidase. *J. Biol. Chem.* 279: 35334-35340.
4. Leary, S.C., et al. 2004. Human SCO1 and SCO2 have independent, cooperative functions in copper delivery to cytochrome c oxidase. *Hum. Mol. Genet.* 13: 1839-1848.
5. Williams, J.C., et al. 2005. Crystal structure of human SCO1: implications for redox signaling by a mitochondrial cytochrome c oxidase assembly protein. *J. Biol. Chem.* 280: 15202-15211.
6. Horng, Y.C., et al. 2005. Human SCO1 and SCO2 function as copper-binding proteins. *J. Biol. Chem.* 280: 34113-34122.
7. Horvath, R., et al. 2005. Congenital cataract, muscular hypotonia, developmental delay and sensori-neural hearing loss associated with a defect in copper metabolism. *J. Inherit. Metab. Dis.* 28: 479-492.

CHROMOSOMAL LOCATION

Genetic locus: SCO1 (human) mapping to 17p13.1.

PRODUCT

SCO1 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see SCO1 shRNA Plasmid (h): sc-61505-SH and SCO1 shRNA (h) Lentiviral Particles: sc-61505-V as alternate gene silencing products.

For independent verification of SCO1 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-61505A, sc-61505B and sc-61505C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

SCO1 siRNA (h) is recommended for the inhibition of SCO1 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

SCO1 (D-1): sc-398001 is recommended as a control antibody for monitoring of SCO1 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor SCO1 gene expression knockdown using RT-PCR Primer: SCO1 (h)-PR: sc-61505-PR (20 μ l, 405 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.