

KIF16B siRNA (m): sc-61592

BACKGROUND

The Kinesins constitute a large family of microtubule-dependent motor proteins which are responsible for the distribution of numerous organelles, vesicles and macromolecular complexes throughout the cell. Individual Kinesin members play crucial roles in cell division, intracellular transport and membrane trafficking events, including endocytosis and transcytosis. KIF16B (kinesin family member 16B), also known as SNX23 (sorting nexin-23) or C20orf23, is a 1,317 amino acid protein that contains one FHA domain, one kinesin-motor domain and one PX domain and belongs to the Kinesin-like protein family. Present in early endosomes at the end of microtubules, KIF16B is thought to function as a microtubule-dependent motor protein that may regulate the motility of early endosomes, thereby mediating the balance between endosomal recycling and degradation. Five isoforms of KIF16B are expressed due to alternative splicing events.

REFERENCES

1. Teasdale, R.D., Loci, D., Houghton, F., Karlsson, L. and Gleeson, P.A. 2001. A large family of endosome-localized proteins related to sorting nexin 1. *Biochem. J.* 358: 7-16.
2. Worby, C.A. and Dixon, J.E. 2002. Sorting out the cellular functions of sorting nexins. *Nat. Rev. Mol. Cell Biol.* 3: 919-931.
3. Hoepfner, S., Severin, F., Cabezas, A., Habermann, B., Runge, A., Gillooly, D., Stenmark, H. and Zerial, M. 2005. Modulation of receptor recycling and degradation by the endosomal kinesin KIF16B. *Cell* 121: 437-450.
4. Miki, H., Okada, Y. and Hirokawa, N. 2005. Analysis of the kinesin superfamily: insights into structure and function. *Trends Cell Biol.* 15: 467-476.
5. Seet, L.F. and Hong, W. 2006. The Phox (PX) domain proteins and membrane traffic. *Biochim. Biophys. Acta* 1761: 878-896.
6. Blatner, N.R., Wilson, M.I., Lei, C., Hong, W., Murray, D., Williams, R.L. and Cho, W. 2007. The structural basis of novel endosome anchoring activity of KIF16B kinesin. *EMBO J.* 26: 3709-3719.
7. Vasilescu, J., Zweitzig, D.R., Denis, N.J., Smith, J.C., Ethier, M., Haines, D.S. and Figeys, D. 2007. The proteomic reactor facilitates the analysis of affinity-purified proteins by mass spectrometry: application for identifying ubiquitinated proteins in human cells. *J. Proteome Res.* 6: 298-305.

CHROMOSOMAL LOCATION

Genetic locus: Kif16b (mouse) mapping to 2 G1.

PRODUCT

KIF16B siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see KIF16B shRNA Plasmid (m): sc-61592-SH and KIF16B shRNA (m) Lentiviral Particles: sc-61592-V as alternate gene silencing products.

For independent verification of KIF16B (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-61592A, sc-61592B and sc-61592C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

KIF16B siRNA (m) is recommended for the inhibition of KIF16B expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

KIF16B (H-6): sc-390309 is recommended as a control antibody for monitoring of KIF16B gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor KIF16B gene expression knockdown using RT-PCR Primer: KIF16B (m)-PR: sc-61592-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.