

Spartin siRNA (h): sc-61601

BACKGROUND

Spartin is a protein that may be involved in microtubule dynamics and endosomal trafficking. The Spartin protein contains a microtubule interacting and trafficking (MIT) molecule domain and is ubiquitously expressed, with highest levels observed in adipose tissue. A frameshift mutation in the Spartin gene (SPG20) causes spastic paraplegia 20, also designated Troyer syndrome, an autosomal recessive form of hereditary spastic paraplegia (HSP). HSP is an inherited neurological disorder characterized by lower extremity weakness and stiffness due to a length-dependent, retrograde axonopathy of corticospinal motor neurons.

REFERENCES

1. Cross, H.E., et al. 1967. The Troyer syndrome. A recessive form of spastic paraplegia with distal muscle wasting. *Arch. Neurol.* 16: 473-485.
2. Neuhauser, G., et al. 1976. Familial spastic paraplegia with distal muscle wasting in the Old Order Amish; atypical Troyer syndrome or new syndrome. *Clin. Genet.* 9: 315-323.
3. Patel, H., et al. 2002. SPG20 is mutated in Troyer syndrome, an hereditary spastic paraplegia. *Nat. Genet.* 31: 347-348.
4. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 275900. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Ciccarelli, F.D., et al. 2003. The identification of a conserved domain in both Spartin and Spastin, mutated in hereditary spastic paraplegia. *Genomics* 81: 437-441.
6. Proukakakis, C., et al. 2004. Troyer syndrome revisited. A clinical and radiological study of a complicated hereditary spastic paraplegia. *J. Neurol.* 251: 1105-1110.
7. Bakowska, J.C., et al. 2005. The Troyer syndrome (SPG20) protein Spartin interacts with Eps15. *Biochem. Biophys. Res. Commun.* 334: 1042-1048.
8. Winner, B., et al. 2006. Thin corpus callosum and amyotrophy in spastic paraplegia-case report and review of literature. *Clin. Neurol. Neurosurg.* 108: 692-698.

CHROMOSOMAL LOCATION

Genetic locus: SPG20 (human) mapping to 13q13.3.

PRODUCT

Spartin siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Spartin shRNA Plasmid (h): sc-61601-SH and Spartin shRNA (h) Lentiviral Particles: sc-61601-V as alternate gene silencing products.

For independent verification of Spartin (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-61601A, sc-61601B and sc-61601C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Spartin siRNA (h) is recommended for the inhibition of Spartin expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

Spartin (D-4): sc-271888 is recommended as a control antibody for monitoring of Spartin gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Spartin gene expression knockdown using RT-PCR Primer: Spartin (h)-PR: sc-61601-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.