



PRSS50 siRNA (h): sc-61730

BACKGROUND

Testis-specific protease 50 (TSP50) is a 385 amino acid protein that is only expressed in the nucleus and cytoplasm of spermatocytes, but not in mature sperm of testis or in malignant mammary epithelial cells. TSP50 is similar in enzymatic structure to many serine proteases, however, in the most critical catalytic triad, the serine has been replaced by a threonine. This unique structure of TSP50 indicates that it may have a specialized substrate. TSP50 is usually methylated, but expression of the abnormally hypomethylated form of TSP50 may lead to breast and ovarian cancer. This suggests that while the normal function of TSP50 is to play a role in spermatogenesis and human regeneration, the upregulation or hypomethylation of TSP50 may be oncogenic.

REFERENCES

1. Yuan, L., Shan, J., De Risi, D., Broome, J., Lovecchio, J., Gal, D., Vinciguerra, V. and Xu, H.P. 1999. Isolation of a novel gene, TSP50, by a hypomethylated DNA fragment in human breast cancer. *Cancer Res.* 59: 3215-3221.
2. Shan, J., Yuan, L., Xiao, Q., Chiorazzi, N., Budman, D., Teichberg, S. and Xu, H.P. 2002. TSP50, a possible protease in human testes, is activated in breast cancer epithelial cells. *Cancer Res.* 62: 290-294.
3. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 607950. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
4. Takano, N., Matsui, H. and Takahashi, T. 2004. TESSP-1: a novel serine protease gene expressed in the spermatogonia and spermatocytes of adult mouse testes. *Mol. Reprod. Dev.* 70: 1-10.
5. Xu, H.P., Yuan, L., Shan, J. and Feng, H. 2004. Localization and expression of TSP50 protein in human and rodent testes. *Urology* 64: 826-832.
6. Looijenga, L.H., Hersmus, R., Gillis, A.J., Pfundt, R., Stoop, H.J., van Gorp, R.J., Veltman, J., Beverloo, H.B., van Drunen, E., van Kessel, A.G., Pera, R.R., Schneider, D.T., Summersgill, B., Shipley, J., et al. 2006. Genomic and expression profiling of human spermatocytic seminomas: primary spermatocyte as tumorigenic precursor and DMRT1 as candidate chromosome 9 gene. *Cancer Res.* 66: 290-302.
7. Wei, J., Liu, Y., Yang, S., Xu, J., Kong, H., Han, B., Bao, Y., Wu, Y., Yin, W., Li, W., Yan, G., Luo, G., Xu, H.P., Li, Y. and Yang, B. 2006. Screening of single-chain variable fragments against TSP50 from a phage display antibody library and their expression as soluble proteins. *J. Biomol. Screen.* 11: 546-552.

CHROMOSOMAL LOCATION

Genetic locus: PRSS50 (human) mapping to 3p21.31.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

PRSS50 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see PRSS50 shRNA Plasmid (h): sc-61730-SH and PRSS50 shRNA (h) Lentiviral Particles: sc-61730-V as alternate gene silencing products.

For independent verification of PRSS50 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-61730A, sc-61730B and sc-61730C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

PRSS50 siRNA (h) is recommended for the inhibition of PRSS50 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor PRSS50 gene expression knockdown using RT-PCR Primer: PRSS50 (h)-PR: sc-61730-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.