

# Vigilin siRNA (h): sc-61788

## BACKGROUND

Vigilin, a K homology (KH) protein, is found in the nucleus and cytoplasm of all eukaryotic species. Vigilin contains a unique structure of 14 to 15 consecutively arranged KH domains, which function to mediate RNA-protein binding. Expression of the gene encoding Vigilin, which maps to chromosome 2q37.3, is essential for cell viability. Vigilin is active in heterochromatin formation and cytoplasmic mRNA decay, and can be a useful marker for translational activity. The 80S ribosome co-localizes with Vigilin, which interacts with the ribosomal complex through its C-terminal domain, suggesting its role in the link between tRNA-export and the channeled tRNA-cycle on ribosomes. Intra-cellular cholesterol upregulates Vigilin expression, and the protein specifically binds to high density lipoprotein molecules to function in the removal of excess cellular cholesterol.

## REFERENCES

- McKnight, G.L., et al. 1992. Cloning and expression of a cellular high density lipoprotein-binding protein that is up-regulated by cholesterol loading of cells. *J. Biol. Chem.* 267: 12131-12141.
- Xia, Y.R., et al. 1993. Localization of the gene for high-density lipoprotein binding protein to human chromosome 2q37. *Genomics* 16: 524-525.
- Plenz, G., et al. 1994. The human vigilin gene: identification, chromosomal localization and expression pattern. *Hum. Genet.* 93: 575-582.
- Kügler, S., et al. 1996. Two additional 5' exons in the human gene and provide the structural basis for differential routes of gene expression. *Eur. J. Biochem.* 238: 410-417.
- Kügler, S., et al. 1996. Vigilin contains a functional nuclear localisation sequence and is present in both the cytoplasm and the nucleus. *FEBS Lett.* 382: 330-334.
- Dodson, R.E., et al. 1997. Vigilin, a ubiquitous protein with 14 K homology domains, is the estrogen-inducible vitellogenin mRNA 3'-untranslated region-binding protein. *J. Biol. Chem.* 272: 12249-12252.
- Chiu, D.S., et al. 1998. High-density lipoprotein-binding protein (HBPI)/Vigilin is expressed in human atherosclerotic lesions and colocalizes with apolipoprotein E. *Arterioscler. Thromb. Vasc. Biol.* 17: 2350-2358.

## CHROMOSOMAL LOCATION

Genetic locus: HDLBP (human) mapping to 2q37.3.

## PRODUCT

Vigilin siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Vigilin shRNA Plasmid (h): sc-61788-SH and Vigilin shRNA (h) Lentiviral Particles: sc-61788-V as alternate gene silencing products.

For independent verification of Vigilin (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-61788A, sc-61788B and sc-61788C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

Vigilin siRNA (h) is recommended for the inhibition of Vigilin expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

Vigilin (H-3): sc-271523 is recommended as a control antibody for monitoring of Vigilin gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Vigilin gene expression knockdown using RT-PCR Primer: Vigilin (h)-PR: sc-61788-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.