

# Na<sup>+</sup>/K<sup>+</sup>-ATPase β3 siRNA (h): sc-62002

## BACKGROUND

The ubiquitously expressed sodium/potassium-ATPase (Na<sup>+</sup>/K<sup>+</sup>-ATPase) exists as a oligomeric plasma membrane complex that couples the hydrolysis of one molecule of ATP to the importation of three Na<sup>+</sup> ions and two K<sup>+</sup> ions against their respective electrochemical gradients. As a member of the P-type family of ion motives, Na<sup>+</sup>/K<sup>+</sup>-ATPase plays a critical role in maintaining cellular volume, resting membrane potential and Na<sup>+</sup>-coupled solute transport. Multiple isoforms of three subunits, α, β and γ, comprise the Na<sup>+</sup>/K<sup>+</sup>-ATPase oligomer. The α subunit contains the binding sites for ATP and the cations; the glycosylated β subunit ensures correct folding and membrane insertion of the α subunits. The small γ subunit co-localizes with the α subunit in nephron segments, where it increases the affinity of Na<sup>+</sup>/K<sup>+</sup>-ATPase for ATP. The β subunit, but not the γ subunit, is essential for normal activity of Na<sup>+</sup>/K<sup>+</sup>-ATPase.

## REFERENCES

1. Hardwicke, P.M., et al. 1981. A proteolipid associated with Na<sup>+</sup>/K<sup>+</sup>-ATPase is not essential for ATPase activity. *Biochem. Biophys. Res. Commun.* 102: 250-257.
2. Ackermann, U., et al. 1990. Mutual dependence of Na<sup>+</sup>/K<sup>+</sup>-ATPase α- and β-subunits for correct post-translational processing and intracellular transport. *FEBS Lett.* 269: 105-108.
3. McDonough, A.A., et al. 1990. The sodium pump needs its β subunit. *FASEB J.* 4: 1598-1605.
4. Pedemonte, C.H., et al. 1990. Chemical modification as an approach to elucidation of sodium pump structure-function relations. *Am. J. Physiol.* 258: C1-C23.
5. Mercer, R.W., et al. 1993. Molecular cloning and immunological characterization of the γ-polypeptide, a small protein associated with Na,K-ATPase. *J. Cell Biol.* 121: 579-586.
6. DeTomaso, A.W., et al. 1993. Expression, targeting, and assembly of functional Na,K-ATPase polypeptides in baculovirus-infected insect cells. *J. Biol. Chem.* 268: 1470-1478.
7. Lingrel, J.B., et al. 1994. Na<sup>+</sup>/K<sup>+</sup>-ATPase. *J. Biol. Chem.* 269: 19659-19662.

## CHROMOSOMAL LOCATION

Genetic locus: ATP1B3 (human) mapping to 3q23.

## PRODUCT

Na<sup>+</sup>/K<sup>+</sup>-ATPase β3 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μM solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Na<sup>+</sup>/K<sup>+</sup>-ATPase β3 shRNA Plasmid (h): sc-62002-SH and Na<sup>+</sup>/K<sup>+</sup>-ATPase β3 shRNA (h) Lentiviral Particles: sc-62002-V as alternate gene silencing products.

For independent verification of Na<sup>+</sup>/K<sup>+</sup>-ATPase β3 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-62002A, sc-62002B and sc-62002C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μl of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μl of RNase-free water makes a 10 μM solution in a 10 μM Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

Na<sup>+</sup>/K<sup>+</sup>-ATPase β3 siRNA (h) is recommended for the inhibition of Na<sup>+</sup>/K<sup>+</sup>-ATPase β3 expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μM in 66 μl. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

Na<sup>+</sup>/K<sup>+</sup>-ATPase β3 (46): sc-135998 is recommended as a control antibody for monitoring of Na<sup>+</sup>/K<sup>+</sup>-ATPase β3 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgGκ BP-FITC: sc-516140 or m-IgGκ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Na<sup>+</sup>/K<sup>+</sup>-ATPase β3 gene expression knockdown using RT-PCR Primer: Na<sup>+</sup>/K<sup>+</sup>-ATPase β3 (h)-PR: sc-62002-PR (20 μl). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.