



Krox-26 siRNA (h): sc-62530

BACKGROUND

Krox-26, also known as ZNF22 (zinc finger protein 22), HKR-T1 or ZNF422, is a 224 amino acid protein that localizes to the nucleus and contains five CDH2-type zinc fingers. One of several members of the Krüppel C₂H₂-type zinc-finger protein family, Krox-26 binds to DNA and is thought to be involved in transcriptional regulation events related to tooth formation, both in embryonic and adult tissue. The gene encoding Krox-26 maps to human chromosome 10, which houses over 1,200 genes and comprises nearly 4.5% of the human genome. Defects in some of the genes that map to chromosome 10 are associated with Charcot-Marie Tooth disease, Jackson-Weiss syndrome, Usher syndrome, nonsyndromic deafness, Wolman's syndrome, Cowden syndrome, multiple endocrine neoplasia type 2 and porphyria.

REFERENCES

- Thiesen, H.J. 1990. Multiple genes encoding zinc finger domains are expressed in human T cells. *New Biol.* 2: 363-374.
- Bray, P., et al. 1991. Characterization and mapping of human genes encoding zinc finger proteins. *Proc. Natl. Acad. Sci. USA* 88: 9563-9567.
- Wu, B.Y., et al. 1992. Isolation of a cDNA clone encoding a zinc finger protein highly expressed in T-leukemia lines. *Blood* 80: 2571-2576.
- Rousseau-Merck, M.F., et al. 1992. A cluster of expressed zinc finger protein genes in the pericentromeric region of human chromosome 10. *Genomics* 13: 845-848.
- Lichter, P., et al. 1992. Clustering of C₂-H₂ zinc finger motif sequences within telomeric and fragile site regions of human chromosomes. *Genomics* 13: 999-1007.
- Love, D.R., et al. 1993. A polymorphic dinucleotide repeat at the ZNF22 locus. *Hum. Mol. Genet.* 2: 491.
- Online Mendelian Inheritance in Man, OMIM™. 1998. Johns Hopkins University, Baltimore, MD. MIM Number: 194529. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Gao, Y., et al. 2003. The human KROX-26/ZNF22 gene is expressed at sites of tooth formation and maps to the locus for permanent tooth agenesis (He-Zhao deficiency). *J. Dent. Res.* 82: 1002-1007.

CHROMOSOMAL LOCATION

Genetic locus: ZNF22 (human) mapping to 10q11.21.

PRODUCT

Krox-26 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Krox-26 shRNA Plasmid (h): sc-62530-SH and Krox-26 shRNA (h) Lentiviral Particles: sc-62530-V as alternate gene silencing products.

For independent verification of Krox-26 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-62530A, sc-62530B and sc-62530C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Krox-26 siRNA (h) is recommended for the inhibition of Krox-26 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Krox-26 gene expression knockdown using RT-PCR Primer: Krox-26 (h)-PR: sc-62530-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.