

LASS5 siRNA (h): sc-62551

BACKGROUND

The LASS (longevity assurance homolog) family members are highly conserved from yeasts to mammals. Six members of this family of proteins have been characterized (LASS1, LASS2, LASS3, LASS4, LASS5 and LASS6) and they are all involved in sphingolipid synthesis. LASS5, also called TRH4, is a 392 amino acid endoplasmic reticulum, multi-pass membrane protein. On the luminal side of the endoplasmic reticulum membrane, the N-terminal asparagine residue is glycosylated. Functioning as a bona fide (dihydro)ceramide synthase, LASS5 increases the levels of short ceramide species, such as C14:0- and C16:0-ceramides as well as those selectively enriched in palmitic acid. LASS5 is the most abundantly expressed longevity assurance homolog in lung epithelia where it functions as a regulator of PtdCho metabolism. In cells deficient for CLN9, as observed in neuronal ceroid lipofuscinosis (NCL) or Batten disease, LASS5 can increase ceramide levels and partially correct growth and apoptosis.

REFERENCES

1. Riebeling, C., Allegood, J.C., Wang, E. and Futerman, A.H. 2003. Two mammalian longevity assurance gene (LAG1) family members, trh1 and trh, regulate dihydroceramide synthesis using different fatty acyl-CoA donors. *J. Biol. Chem.* 278: 43452-43459.
2. Spassieva, S., Seo, J.G., Jiang, J.C., Bielawski, J., Alvarez-Vasquez, F., Jazwinski, S.M., Hannun, Y.A. and Obeid, L.M. 2006. Necessary role for the Lag1p motif in (dihydro)ceramide synthase activity. *J. Biol. Chem.* 281: 33931-33938.
3. Mizutani, Y., Kihara, A. and Igarashi, Y. 2005. Mammalian LASS6 and its related family members regulate synthesis of specific ceramides. *Biochem. J.* 390: 263-271.
4. Xu, Z., Zhou, J., McCoy, D.M. and Mallampalli, R.K. 2005. LASS5 is the predominant ceramide synthase isoform involved in *de novo* sphingolipid synthesis in lung epithelia. *J. Lipid Res.* 46: 1229-1238.
5. Lahiri, S. and Futerman, A.H. 2005. LASS5 is a bona fide dihydroceramide synthase that selectively utilizes palmitoyl-CoA as acyl donor. *J. Biol. Chem.* 280: 33735-33738.
6. Schulz, A., Mousallem, T., Venkataramani, M., Persaud-Sawin, D.A., Zucker, A., Luberto, C., Bielawska, A., Bielawski, J., Holthuis, J.C., Jazwinski, S.M., Kozhaya, L., Dbaiho, G.S. and Boustany, R.M. 2006. The CLN9 protein, a regulator of dihydroceramide synthase. *J. Biol. Chem.* 281: 2784-2794.

CHROMOSOMAL LOCATION

Genetic locus: CERS5 (human) mapping to 12q13.12.

PRODUCT

LASS5 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see LASS5 shRNA Plasmid (h): sc-62551-SH and LASS5 shRNA (h) Lentiviral Particles: sc-62551-V as alternate gene silencing products.

For independent verification of LASS5 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-62551A, sc-62551B and sc-62551C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

LASS5 siRNA (h) is recommended for the inhibition of LASS5 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor LASS5 gene expression knockdown using RT-PCR Primer: LASS5 (h)-PR: sc-62551-PR (20 μ l, 474 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.