

MTHFS siRNA (h): sc-62649

BACKGROUND

MTHFS (5-formyltetrahydrofolate cyclo-ligase) is a cytosolic protein involved in the formate metabolic process. MTHFS can account for up to thirty percent of a cell's cytoplasmic folate derivatives. MTHFS, with a magnesium cofactor, catalyzes the ATP-dependent reaction that reduces 5-formyltetrahydrofolate to 5,10-methylenetetrahydrofolate. Folate is a necessary molecule for DNA replication and a deficiency in folate can lead to numerous disease states. During DNA replication 5,10-methylenetetrahydrofolate is oxidized and MTHFS is responsible for resetting the molecule.

REFERENCES

1. Anguera, M.C., Liu, X. and Stover, P.J. 2004. Cloning, expression, and purification of 5,10-methylenetetrahydrofolate synthetase from *Mus musculus*. *Protein Expr. Purif.* 35: 276-283.
2. Chen, S., Yakunin, A.F., Proudfoot, M., Kim, R. and Kim, S.H. 2005. Structural and functional characterization of a 5,10-methylenetetrahydrofolate synthetase from *Mycoplasma pneumoniae* (GI: 13508087). *Proteins* 61: 433-443.
3. Anguera, M.C. and Stover, P.J. 2006. Methylenetetrahydrofolate synthetase is a high-affinity catecholamine-binding protein. *Arch. Biochem. Biophys.* 455: 175-187.
4. Field, M.S., Szebenyi, D.M. and Stover, P.J. 2006. Regulation of *de novo* purine biosynthesis by methylenetetrahydrofolate synthetase in neuroblastoma. *J. Biol. Chem.* 281: 4215-4221.
5. Lee, K.M., Lan, Q., Krickler, A., Purdue, M.P., Grulich, A.E., Vajdic, C.M., Turner, J., Whitby, D., Kang, D., Chanock, S., Rothman, N. and Armstrong, B.K. 2007. One-carbon metabolism gene polymorphisms and risk of non-Hodgkin lymphoma in Australia. *Hum. Genet.* 122: 525-533.
6. Matakidou, A., El Gatta, R., Rudd, M.F., Webb, E.L., Bridle, H., Eisen, T. and Houlston, R.S. 2007. Prognostic significance of folate metabolism polymorphisms for lung cancer. *Br. J. Cancer* 97: 247-252.

CHROMOSOMAL LOCATION

Genetic locus: MTHFS (human) mapping to 15q25.1.

PRODUCT

MTHFS siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see MTHFS shRNA Plasmid (h): sc-62649-SH and MTHFS shRNA (h) Lentiviral Particles: sc-62649-V as alternate gene silencing products.

For independent verification of MTHFS (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-62649A, sc-62649B and sc-62649C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

MTHFS siRNA (h) is recommended for the inhibition of MTHFS expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

MTHFS (G-8): sc-398056 is recommended as a control antibody for monitoring of MTHFS gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor MTHFS gene expression knockdown using RT-PCR Primer: MTHFS (h)-PR: sc-62649-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.