

StARD7 siRNA (h): sc-63078

BACKGROUND

The StARD (steroidogenic acute regulatory protein-related lipid transfer (START) domain containing) family of proteins is comprised of 15 different members. All members contain the characteristic START domain and are believed to play key roles in the metabolism and transport of lipids. The StARD proteins are grouped into six subfamilies based on their START domain sequences. PC-TP (StARD2), StARD7, StARD10 and GPBP (StARD11) constitute one subfamily, namely the STARD2/PCTP group. StARD7, also known as GTT1 (gestational trophoblastic tumor 1), is a widely expressed protein with a possible function in phospholipid transport. The sequence of its START domain suggests that StARD7 may bind to charged lipids. StARD7 is overexpressed in a number of cancer cell lines and thus may be involved in tumor signaling mediated by phospholipids.

REFERENCES

1. Soccio, R.E. and Breslow, J.L. 2003. StAR-related lipid transfer (START) proteins: mediators of intracellular lipid metabolism. *J. Biol. Chem.* 278: 22183-22186.
2. Durand, S., Angeletti, S. and Genti-Raimondi, S. 2004. GTT1/StarD7, a novel phosphatidylcholine transfer protein-like highly expressed in gestational trophoblastic tumour: cloning and characterization. *Placenta* 25: 37-44.
3. Angeletti, S., Maggio, B. and Genti-Raimondi, S. 2004. Surface activity and interaction of StarD7 with phospholipid monolayers. *Biochem. Biophys. Res. Commun.* 314: 181-185.
4. Soccio, R.E., Adams, R.M., Maxwell, K.N. and Breslow, J.L. 2005. Differential gene regulation of StarD4 and StarD5 cholesterol transfer proteins. Activation of StarD4 by sterol regulatory element-binding protein-2 and StarD5 by endoplasmic reticulum stress. *J. Biol. Chem.* 280: 19410-19418.
5. Alpy, F. and Tomasetto, C. 2005. Give lipids a START: the StAR-related lipid transfer (START) domain in mammals. *J. Cell Sci.* 118: 2791-2801.
6. Olaiyoye, M.A., Vehring, S., Müller, P., Herrmann, A., Schiller, J., Thiele, C., Lindeman, G.J., Visvader, J.E. and Pomorski, T. 2005. StarD10, a START domain protein overexpressed in breast cancer, functions as a phospholipid transfer protein. *J. Biol. Chem.* 280: 27436-27442.
7. Tian, J., Zhang, S., Liu, Z., Zhuang, Y., Wang, Y. and Jiang, S. 2007. Characterization and tissue-specific expression of phosphatidylcholine transfer protein gene from amphioxus *Branchiostoma belcheri*. *Cell Tissue Res.* 330: 53-61.
8. Lee, H.S., Park, M.H., Yang, S.J., Park, K.C., Kim, N.S., Kim, Y.S., Kim, D.I., Yoo, H.S., Choi, E.J. and Yeom, Y.I. 2007. Novel candidate targets of Wnt/ β -catenin signaling in hepatoma cells. *Life Sci.* 80: 690-698.

CHROMOSOMAL LOCATION

Genetic locus: STARD7 (human) mapping to 2q11.2.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

StARD7 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see StARD7 shRNA Plasmid (h): sc-63078-SH and StARD7 shRNA (h) Lentiviral Particles: sc-63078-V as alternate gene silencing products.

For independent verification of StARD7 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-63078A, sc-63078B and sc-63078C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

StARD7 siRNA (h) is recommended for the inhibition of StARD7 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor StARD7 gene expression knockdown using RT-PCR Primer: StARD7 (h)-PR: sc-63078-PR (20 μ l, 509 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.