

CD71 (3H3074): sc-70778

BACKGROUND

CD71, also known as the transferrin receptor (TFR), is a type II membrane glycoprotein that exists as a disulfide-linked homodimer of two identical subunits. CD71 binds to two molecules of transferrin and a serum iron-transport protein and directs the cellular uptake of iron via receptor-mediated endocytosis. CD71 is expressed, typically at high levels, on all proliferating cells, reticulocytes and erythroid precursors. It is not expressed on resting leukocytes, but is upregulated upon activation of lymphocytes, monocytes and macrophages. CD71 is also found on most dividing cells and on brain endothelium. A second transferrin receptor, TFR2, also mediates the uptake of transferrin-bound iron. TFR2 is a homodimer of two subunits and is highly expressed in liver as well as in hepatocytes and erythroid precursors. Mutations in the TFR2 gene result in hereditary hemochromatosis type III (HFE3), an iron overloading disorder predominant in Caucasians.

REFERENCES

1. Lesley, J., et al. 1984. Expression of transferrin receptor on murine hematopoietic progenitors. *Cell. Immunol.* 83: 14-25.
2. McClelland, A., et al. 1984. The human transferrin receptor gene: genomic organization, and the complete primary structure of the receptor deduced from a cDNA sequence. *Cell* 39: 267-274.
2. Lesley, J.F., et al. 1986. Inhibition of cell growth by monoclonal anti-transferrin receptor antibodies. *Mol. Cell. Biol.* 5: 1814-1821.
3. Kemp, J.D., et al. 1987. Role of the transferrin receptor in lymphocyte growth: a rat IgG monoclonal antibody against the murine transferrin receptor produces highly selective inhibition of T and B cell activation protocols. *J. Immunol.* 138: 2422-2426.
4. Sauvage, C.A., et al. 1987. Effects of monoclonal antibodies that block transferrin receptor function on the *in vivo* growth of a syngeneic murine leukemia. *Cancer Res.* 47: 747-753.
5. Kemp, J.D., et al. 1989. Inhibition of lymphocyte activation with anti-transferrin receptor Mabs: a comparison of three reagents and further studies of their range of effects and mechanism of action. *Cell. Immunol.* 122: 218-230.
6. Testa, U., et al. 1993. The transferrin receptor. *Crit. Rev. Oncog.* 4: 241-276.
7. Brekelmans, P., et al. 1995. Transferrin receptor expression as a marker of immature cycling thymocytes in the mouse. *Cell. Immunol.* 159: 331-339.
8. Pasquier, B., et al. 2004. Differential expression and function of IgA receptors (CD89 and CD71) during maturation of dendritic cells. *J. Leukoc. Biol.* 76: 1134-1141.

STORAGE

Store at 4° C, **DO NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: Tfrc (mouse) mapping to 16 B3.

SOURCE

CD71 (3H3074) is a rat monoclonal antibody raised against concanavalin A-stimulated spleen cells of mouse origin.

PRODUCT

Each vial contains 200 µg IgG_{2a} in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

CD71 (3H3074) is recommended for detection of CD71 of mouse and rat origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and flow cytometry (1 µg per 1 x 10⁶ cells).


Suitable for use as control antibody for CD71 siRNA (m): sc-37071, CD71 shRNA Plasmid (m): sc-37071-SH and CD71 shRNA (m) Lentiviral Particles: sc-37071-V.

Molecular Weight of CD71: 85-95 kDa.


Molecular Weight of CD71 dimer: 190 kDa.

Positive Controls: M1 whole cell lysate: sc-364782, MH-S whole cell lysate: sc-364785 or RAW 264.7 whole cell lysate: sc-2211.

DATA



CD71 (3H3074): sc-70778. Western blot analysis of CD71 expression in M1 (A) and mouse PBL (B) whole cell lysates.



CD71 (3H3074): sc-70778. Western blot analysis of CD71 expression in MH-S whole cell lysate.

RESEARCH USE

For research use only, not for use in diagnostic procedures.



See **CD71 (YTA 74.4): sc-59112** for CD71 antibody conjugates, including AC, HRP, FITC, PE, and Alexa Fluor® 488, 546, 594, 647, 680 and 790.