

## CD74 (2X31): sc-70781



The Power to Question

## BACKGROUND

The human histocompatibility leukocyte antigen (HLA) class II-associated invariant chain is composed of at least four polypeptides. One of these polypeptide chains is expressed as a membrane-bound subunit and has been designated CD74. The loading of peptide onto the class II MHC protein (MHC II) appears to be regulated by CD74, which associates with MHC II during its migration to the endosomal compartment, where class II binds peptide. CD74 is expressed by cells of both T lymphocyte and B lymphocyte lineages. In fact, CD74 is broadly expressed in normal B lymphocytes, regardless of their histocompatibility leukocyte antigen (HLA) phenotype, while a subset of peripheral T lymphocytes that are MHC II negative do not express CD74.

## REFERENCES

1. McMaster, W.R., et al. 1979. Identification of Ia glycoproteins in rat thymus and purification from rat spleen. *Eur. J. Immunol.* 9: 426-433.
2. Sarker, A.B., et al. 1992. *Bauhinia purpurea*—a new paraffin section marker for Reed-Sternberg cells of Hodgkin's disease. A comparison with Leu-M1 (CD15), LN2 (CD74), peanut agglutinin and Ber-H2 (CD30). *Am. J. Pathol.* 141: 19-23.

## CHROMOSOMAL LOCATION

Genetic locus: CD74 (human) mapping to 5q32; Cd74 (mouse) mapping to 18 E1.

## SOURCE

CD74 (2X31) is a mouse monoclonal antibody raised against SU-DHL-4 lymphoma cells of human origin.

## PRODUCT

Each vial contains 200 µg IgG<sub>1</sub> kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

## APPLICATIONS

CD74 (2X31) is recommended for detection of CD74 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500), immunohistochemistry (including paraffin-embedded sections) (starting dilution 1:50, dilution range 1:50-1:500) and flow cytometry (1 µg per 1 x 10<sup>6</sup> cells).

Suitable for use as control antibody for CD74 siRNA (h): sc-35023, CD74 siRNA (m): sc-35024, CD74 shRNA Plasmid (h): sc-35023-SH, CD74 shRNA Plasmid (m): sc-35024-SH, CD74 shRNA (h) Lentiviral Particles: sc-35023-V and CD74 shRNA (m) Lentiviral Particles: sc-35024-V.


Molecular Weight of CD74 isoforms: 31-45 kDa.

Positive Controls: BJAB whole cell lysate: sc-2207, Ramos cell lysate: sc-2216 or GA-10 whole cell lysate: sc-364230.


## RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use m-IgGκ BP-FITC: sc-516140 or m-IgGκ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850. 4) Immunohistochemistry: use m-IgGκ BP-HRP: sc-516102 with DAB, 50X: sc-24982 and Immunohistomount: sc-45086, or Organo/Limonene Mount: sc-45087.

## DATA



CD74 (2X31): sc-70781. Western blot analysis of CD74 expression in Ramos (A), GA-10 (B) and BJAB (C) whole cell lysates and human tonsil tissue extract (D).



CD74 (2X31): sc-70781. Immunoperoxidase staining of formalin fixed, paraffin-embedded human appendix tissue showing membrane staining of lymphoid cells (A). Immunoperoxidase staining of formalin fixed, paraffin-embedded human lymph node tissue showing membrane staining of cells in germinal center and cells in non-germinal center (B).

## SELECT PRODUCT CITATIONS

1. Gai, J.W., et al. 2018. Expression of CD74 in bladder cancer and its suppression in association with cancer proliferation, invasion and angiogenesis in HT-1376 cells. *Oncol. Lett.* 15: 7631-7638.

## STORAGE

Store at 4° C, \*\*DO NOT FREEZE\*\*. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.



See **CD74 (LN-2): sc-6262** for CD74 antibody conjugates, including AC, HRP, FITC, PE, and Alexa Fluor® 488, 546, 594, 647, 680 and 790.