



## RP105 (3H1961): sc-71967

### BACKGROUND

RP105 is a mouse B cell surface molecule that transmits a growth-promoting signal and is implicated in the life/death decision of B cells. RP105 has tandem repeats of a leucine-rich motif in its extracellular domain that are presumed to be involved in protein-protein interactions. The amino acid sequence of human RP105 is highly homologous to that of mouse RP105; human RP105 shares 74% identity with the mouse protein, as well as the leucine-rich motif. Surface expression of RP105 is enhanced in the presence of MD-1, although this expression is restricted to CD19-positive B cells. RP105 demonstrates predominant expression on mature B cells in mantle zones; very little expression is observed in germinal centers.

### REFERENCES

1. Gay, N.J. and Keith, F.J. 1991. *Drosophila* toll and IL-1 receptor. Nature 351: 355-356.
2. Miura, Y., Miyake, K., Yamashita, Y., Shimazu, R., Copeland, N.G., Gilbert, D.J., Jenkins, N.A., Inazawa, J., Abe, T. and Kimoto, M. 1996. Molecular cloning of a human RP105 homologue and chromosomal localization of the mouse and human RP105 genes (Ly64 and LY64). Genomics 38: 299-304.
3. Rock, F.L., Hardiman, G., Timans, J.C., Kastelein, R.A. and Bazan, J.F. 1998. A family of human receptors structurally related to *Drosophila* toll. Proc. Natl. Acad. Sci. USA 95: 588-593.
4. Miura, Y., Shimazu, R., Miyake, K., Akashi, S., Ogata, H., Yamashita, Y., Narisawa, Y. and Kimoto, M. 1998. RP105 is associated with MD-1 and transmits an activation signal in human B cells. Blood 92: 2815-2822.
5. Miyake, K., Shimazu, R., Kondo, J., Niki, T., Akashi, S., Ogata, H., Yamashita, Y., Miura, Y. and Kimoto, M. 1998. Mouse MD-1, a molecule that is physically associated with RP105 and positively regulates its expression. J. Immunol. 161: 1348-1353.
6. Brightbill, H.D., Libraty, D.H., Krutzik, S.R., Yang, R.B., Belisle, J.T., Bleharski, J.R., Maitland, M., Norgard, M.V., Plevy, S.E., Smale, S.T., Brennan, P.J., Bloom, B.R., Godowski, P.J. and Modlin, R.L. 1999. Host defense mechanisms triggered by microbial lipoproteins through toll-like receptors. Science 285: 732-736.
7. Medzhitov, R., Preston-Hurlburt, P. and Janeway, C.A., Jr. 2000. A human homologue of the *Drosophila* toll protein signals activation of adaptive immunity. Nature 388: 394-397.
8. Chuang, T.H. and Ulevitch, R.J. 2000. Cloning and characterization of a subfamily of human toll-like receptors: hTLR7, hTLR8, hTLR9. Eur. Cytokine Netw. 11: 372-378.
9. Miyake, K., Ogata, H., Nagai, Y., Akashi, S. and Kimoto, M. 2000. Innate recognition of lipopolysaccharide by toll-like receptor 4/MD-2 and RP105/MD-1. J. Endotoxin Res. 6: 389-391.

### CHROMOSOMAL LOCATION

Genetic locus: Cd180 (mouse) mapping to 13 D1.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.

### SOURCE

RP105 (3H1961) is a rat monoclonal antibody raised against transfected cells expressing full length RP105 (also known as CD180) of mouse origin.

### PRODUCT

Each vial contains 200 µg IgG<sub>2a</sub> in 1.0 ml PBS with < 0.1% sodium azide and 0.1% gelatin.

### APPLICATIONS

RP105 (3H1961) is recommended for detection of RP105 of mouse origin by immunoprecipitation [1-2 µg per 100-500 µg of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and flow cytometry (1 µg per 1 x 10<sup>6</sup> cells).

Suitable for use as control antibody for RP105 siRNA (m): sc-40253, RP105 shRNA Plasmid (m): sc-40253-SH and RP105 shRNA (m) Lentiviral Particles: sc-40253-V.

Molecular Weight of RP105: 95-105 kDa.

### STORAGE

Store at 4° C, **\*\*DO NOT FREEZE\*\***. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.