



RT1-B α (6A603): sc-71979

BACKGROUND

RT1 is the the major histocompatibility complex (MHC) in the rat which contains genes that code for two class II histocompatibility antigens. RT1-B is an antigen of the RT1 complex. It is a protein dimer consisting of an α and β glycoprotein chain and is homologous to I-A and I-E genes, respectively, in the H-2 complex of the mouse. MHC Class II antigens are useful in studying T helper cell interaction with class II positive antigen presenting cells (dendritic cells, B cells, macrophages) and offer new possibilities for studying the development of T helper cells since these antibodies also stain stromal cells in the thymus.

REFERENCES

1. Holowachuk, E.W., Greer, M.K. and Martin, D.R. 1988. Elevated mRNA levels of major histocompatibility complex class II genes in lymphocytes of autoimmune BB rats. *Diabetes* 37: 1637-1640.
2. Mizuno, K., Inomata, T., Tsuchimoto, S., Matsuno, Y., Niiyama, T., Fujii, H., Natori, T. and Aizawa, M. 1988. Ir gene for bovine insulin in the rat maps to RT1-B β . *J. Immunogenet.* 14: 159-162.
3. Fujii, H., Kakinuma, M., Yoshiki, T. and Natori, T. 1991. Polymorphism of the class II gene of rat major histocompatibility complex, RT1: partial sequence comparison of the first domain of the RT1-B β 1 alleles. *Immunogenetics* 33: 399-403.
4. Syha-Jedelhauser, J., Wendling, U. and Reske, K. 1991. Complete coding nucleotide sequence of cDNA for the class II RT1-B β I chain of the Lewis rat. *Biochim. Biophys. Acta* 1089: 414-416.
5. Easterfield, A.J., Bradley, J.A. and Bolton, E.M. 2003. Complementary DNA sequences encoding the rat MHC class II RT1-Bu and RT1-Du α and β chains. *Immunogenetics* 55: 344-350.
6. Ettinger, R.A., Moustakas, A.K. and Lobaton, S.D. 2004. Open reading frame sequencing and structure-based alignment of polypeptides encoded by RT1-B β , RT1-B α , RT1-D β and RT1-D α alleles. *Immunogenetics* 56: 585-596.
7. Neumann, J. 2005. Novel antibody tags from the rat lysosomal protein RT1.DM for immunodetection of recombinant proteins. *J. Immunol. Methods* 301: 66-76.
8. Roos, C. and Walter, L. 2005. Considerable haplotypic diversity in the RT1-CE class I gene region of the rat major histocompatibility complex. *Immunogenetics* 56: 773-777.

SOURCE

RT1-B α (6A603) is a mouse monoclonal antibody raised against bone marrow cells of rat origin.

PRODUCT

Each vial contains 100 μ g IgG₁ in 1.0 ml PBS with < 0.1% sodium azide and 0.1% gelatin.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

RT1-B α (6A603) is recommended for detection of RT1-B α of mouse and rat origin by immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

Molecular Weight of RT1-B α : 29 kDa.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.