

nm23-H2 siRNA (r): sc-72195

BACKGROUND

The nm23 gene, a potential suppressor of metastasis, was originally identified by differential hybridization between two murine melanoma sub-lines, one with a high and the second with a low metastatic capacity. Highly metastatic sub-lines exhibit much lower levels of nm23 than less metastatic cells. Based on sequence analysis, nm23 appears highly related to nucleotide diphosphate kinases (NDP). In humans, NDP kinases A and B are identical to two isoforms of human nm23 homologs, namely nm23-H1 and H2, respectively. nm23-H2 is identical in sequence to PuF, a transcription factor that binds to nuclease-hypersensitive elements at positions 142 to 115 of the human C-Myc promoter.

REFERENCES

1. Steeg, P.S., Bevilacqua, G., Kopper, L., Thorgeirsson, U.P., Talmadge, J.E., Liotta, L.A. and Sobel, M.E. 1988. Evidence for a novel gene associated with low tumor metastatic potential. *J. Natl. Cancer Inst.* 80: 200-209.
2. Lacombe, M., Wallet, V., Troll, H. and Véron, M. 1990. Functional cloning of a nucleoside diphosphate kinase from *Dictyostelium discoideum*. *J. Biol. Chem.* 265: 10012-10018.
3. Kimura, N., Shimada, N., Nomura, K. and Watanabe, K. 1990. Isolation and characterization of a cDNA clone encoding rat nucleoside diphosphate kinase. *J. Biol. Chem.* 265: 15744-15749.
4. Stahl, J.A., Leone, A., Rosengard, A.M., Porter, L., King, C.R. and Steeg, P.S. 1991. Identification of a second human nm23 gene, nm23-H2. *Cancer Res.* 51: 445-449.
5. Urano, T., Takamiya, K., Furukawa, K. and Shiki, H. 1992. Molecular cloning and functional expression of the second mouse nm23/NDP kinase gene, nm23-M2. *FEBS Lett.* 309: 358-362.
6. Urano, T., Furukawa, K. and Shiku, H. 1993. Expression of nm23/NDP kinase proteins on the cell surface. *Oncogene* 8: 1371-1376.

CHROMOSOMAL LOCATION

Genetic locus: Nme2 (rat) mapping to 10q26.

PRODUCT

nm23-H2 siRNA (r) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see nm23-H2 shRNA Plasmid (r): sc-72195-SH and nm23-H2 shRNA (r) Lentiviral Particles: sc-72195-V as alternate gene silencing products.

For independent verification of nm23-H2 (r) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-72195A, sc-72195B and sc-72195C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

nm23-H2 siRNA (r) is recommended for the inhibition of nm23-H2 expression in rat cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

nm23-H1/2/3 (E-11): sc-166937 is recommended as a control antibody for monitoring of nm23-H2 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor nm23-H2 gene expression knockdown using RT-PCR Primer: nm23-H2 (r)-PR: sc-72195-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.