



Enterokinase siRNA (m): sc-72234

BACKGROUND

Enterokinase, also known as enteropeptidase or serine protease 7, belongs to the peptidase S1 family and localizes to the intestinal brush border in the proximal small intestine. It exists as a heterodimer of a catalytic light chain (LC) and a non-catalytic heavy chain (HC) linked together by a disulfide bond. Enterokinase HC plays a role in macromolecular substrate recognition and specificity. Duodenase is the serine protease responsible for the release and activation of Enterokinase from its inactive precursor. Active Enterokinase recognizes the target sequence, Asp-Asp-Asp-Asp-Lys, and is responsible for catalyzing the conversion of pancreatic trypsinogen to activated trypsin. Activated trypsin then further activates digestive enzymes such as chymotrypsin, carboxypeptidases, elastases and lipases, releasing them from their inactive precursors. Enterokinase is important for proper digestion of proteins. Improper functioning of Enterokinase may result in congenital enteropeptidase deficiency. This recessively inherited disorder leads to severe protein malabsorption and can result in low serum protein, chronic diarrhea and, in infants, a failure to thrive.

REFERENCES

1. Hadorn, B., et al. 1975. Intestinal Enterokinase deficiency. Occurrence in two sibs and age dependency of clinical expression. Arch. Dis. Child. 50: 277-282.
2. Matsushima, M., et al. 1994. Structural characterization of porcine enteropeptidase. J. Biol. Chem. 269: 19976-19982.
3. Kitamoto, Y., et al. 1994. Enterokinase, the initiator of intestinal digestion, is a mosaic protease composed of a distinctive assortment of domains. Proc. Natl. Acad. Sci. USA 91: 7588-7592.
4. Kitamoto, Y., et al. 1995. cDNA sequence and chromosomal localization of human Enterokinase, the proteolytic activator of trypsinogen. Biochemistry 34: 4562-4568.
5. Lu, D., et al. 1998. Bovine proenteropeptidase is activated by trypsin, and the specificity of enteropeptidase depends on the heavy chain. J. Biol. Chem. 272: 31293-31300.
6. Zamolodchikova, T.S., et al. 2000. Activation of recombinant proenteropeptidase by duodenase. FEBS Lett. 466: 295-299.
7. Holzinger, A., et al. 2001. Mutations in the proenteropeptidase gene are the molecular cause of congenital enteropeptidase deficiency. Am. J. Hum. Genet. 70: 20-25.
8. Song, H.W., et al. 2002. Engineered recombinant enteropeptidase catalytic subunit: effect of N-terminal modification. Arch. Biochem. Biophys. 400: 1-6.
9. Likhareva, V.V., et al. 2003. Hydrolysis by enteropeptidase of nonspecific (model) peptide sequences and possible physiological role of this phenomenon. Vopr. Med. Khim. 48: 561-569.

CHROMOSOMAL LOCATION

Genetic locus: Tmprss15 (mouse) mapping to 16 C3.1.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

Enterokinase siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Enterokinase shRNA Plasmid (m): sc-72234-SH and Enterokinase shRNA (m) Lentiviral Particles: sc-72234-V as alternate gene silencing products.

For independent verification of Enterokinase (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-72234A, sc-72234B and sc-72234C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Enterokinase siRNA (m) is recommended for the inhibition of Enterokinase expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Enterokinase gene expression knockdown using RT-PCR Primer: Enterokinase (m)-PR: sc-72234-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.