



GPB5 siRNA (m): sc-72243

BACKGROUND

Glycoprotein hormone β -5 subunit (GPB5) is a cystine knot-forming polypeptide found primarily in the brain and pituitary with low levels expressed in the retina, testis and skin. GPB5 glycoprotein can form a heterodimer with the α subunit, GPA2, both of which are members of the glycoprotein hormone family. Together with GPA2, GPB5 functions to stimulate the thyroid by activating the THSR (thyroid-stimulating hormone receptor), thereby increasing cAMP production. Overexpression of GPB5 is thought to cause reductions in body weight, proptosis and an elevated T4 cell count.

REFERENCES

1. Nakabayashi, K., Matsumi, H., Bhalla, A., Bae, J., Mosselman, S., Hsu, S.Y. and Hsueh, A.J. 2002. Thyrostimulin, a heterodimer of two new human glycoprotein hormone subunits, activates the thyroid-stimulating hormone receptor. *J. Clin. Invest.* 109: 1445-1452.
2. Hsu, S.Y., Nakabayashi, K., Bhalla, A. 2002. Evolution of glycoprotein hormone subunit genes in bilateral metazoa: identification of two novel human glycoprotein hormone subunit family genes, GPA2 and GPB5. *Mol. Endocrinol.* 16: 1538-1551.
3. Macdonald, L.E., Wortley, K.E., Gowen, L.C., Anderson, K.D., Murray, J.D., Poueymirou, W.T., Simmons, M.V., Barber, D., Valenzuela, D.M., Economides, A.N., Wiegand, S.J., Yancopoulos, G.D., Sleeman, M.W. and Murphy, A.J. 2005. Resistance to diet-induced obesity in mice globally overexpressing OGH/GPB5. *Proc. Natl. Acad. Sci. USA* 102: 2496-2501.
4. Park, J.I., Semyonov, J., Chang, C.L. and Hsu, S.Y. 2005. Conservation of the heterodimeric glycoprotein hormone subunit family proteins and the LGR signaling system from nematodes to humans. *Endocrine* 26: 267-276.
5. Sudo, S., Kuwabara, Y., Park, J.I., Hsu, S.Y. and Hsueh, A.J. 2005. Heterodimeric fly glycoprotein hormone α -2 (GPA2) and glycoprotein hormone β -5 (GPB5) activate fly leucine-rich repeat-containing G protein-coupled receptor-1 (DLGR1) and stimulation of human thyrotropin. *Endocrinology* 146: 3596-3604.

CHROMOSOMAL LOCATION

Genetic locus: Gphb5 (mouse) mapping to 12 C3.

PRODUCT

GPB5 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see GPB5 shRNA Plasmid (m): sc-72243-SH and GPB5 shRNA (m) Lentiviral Particles: sc-72243-V as alternate gene silencing products.

For independent verification of GPB5 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-72243A, sc-72243B and sc-72243C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

GPB5 siRNA (m) is recommended for the inhibition of GPB5 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor GPB5 gene expression knockdown using RT-PCR Primer: GPB5 (m)-PR: sc-72243-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.