

Antiquitin siRNA (m): sc-72508

BACKGROUND

Aldehyde dehydrogenases (ALDHs) mediate the NADP⁺-dependent oxidation of aldehydes into acids and play an important role in the detoxification of alcohol-derived acetaldehyde, as well as in lipid peroxidation and in the metabolism of corticosteroids, biogenic amines and neurotransmitters. Antiquitin, also known as ALDH7A1 (aldehyde dehydrogenase 7 family, member A1), ATQ1, EPD or PDE, is a 539 amino acid protein that is highly expressed in ovary, heart, eye, kidney and ear tissue and plays an important role in the detoxification of lipid peroxidation- and alcohol metabolism-related aldehydes. Mutations in the gene encoding Antiquitin are the cause of pyridoxine-dependent epilepsy (PDE), a rare disorder that is characterized by seizures that begin at infancy and involve muscle rigidity, convulsions and loss of consciousness. Additionally, PDE is associated with poor muscle tone, hypothermia and irritability.

REFERENCES

1. Skvorak, A.B., et al. 1997. An ancient conserved gene expressed in the human inner ear: identification, expression analysis, and chromosomal mapping of human and mouse Antiquitin (ATQ1). *Genomics* 46: 191-199.
2. Mills, P.B., et al. 2006. Mutations in Antiquitin in individuals with pyridoxine-dependent seizures. *Nat. Med.* 12: 307-309.
3. Salomons, G.S., et al. 2007. An intriguing "silent" mutation and a founder effect in Antiquitin (ALDH7A1). *Ann. Neurol.* 62: 414-418.
4. Bok, L.A., et al. 2007. Pyridoxine-dependent seizures in Dutch patients: diagnosis by elevated urinary α -amino adipic semialdehyde levels. *Arch. Dis. Child.* 92: 687-689.
5. Plecko, B., et al. 2007. Biochemical and molecular characterization of 18 patients with pyridoxine-dependent epilepsy and mutations of the Antiquitin (ALDH7A1) gene. *Hum. Mutat.* 28: 19-26.
6. Kanno, J., et al. 2007. Allelic and non-allelic heterogeneities in pyridoxine dependent seizures revealed by ALDH7A1 mutational analysis. *Mol. Genet. Metab.* 91: 384-389.

CHROMOSOMAL LOCATION

Genetic locus: Aldh7a1 (mouse) mapping to 18 D3.

PRODUCT

Antiquitin siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Antiquitin shRNA Plasmid (m): sc-72508-SH and Antiquitin shRNA (m) Lentiviral Particles: sc-72508-V as alternate gene silencing products.

For independent verification of Antiquitin (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-72508A, sc-72508B and sc-72508C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Antiquitin siRNA (m) is recommended for the inhibition of Antiquitin expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

Antiquitin (D-7): sc-376906 is recommended as a control antibody for monitoring of Antiquitin gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Antiquitin gene expression knockdown using RT-PCR Primer: Antiquitin (m)-PR: sc-72508-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.