

## FKBP11 siRNA (h): sc-75021

### BACKGROUND

The immunophilins are a highly conserved family of *cis-trans* peptidyl-prolyl isomerases that bind to and mediate the effects of immunosuppressive drugs, such as cyclosporin, FK506 and rapamycin. Immunophilins have also been implicated in protein folding and trafficking within the endoplasmic reticulum (ER). FKBP11 (FK506-binding protein 11), also known as FKBP19 or peptidyl-prolyl *cis-trans* isomerase FKBP11, is a 201 amino acid single-pass membrane protein belonging to the FKBP-type PPlase family, a group of proteins known to catalyze the folding of proline-containing polypeptides. Containing one PPlase FKBP-type domain, FKBP11 is expressed in secretory tissues such as pancreas, pituitary, stomach, lymph node and salivary gland, and is encoded by a gene that maps to human chromosome 12q13.12. FK506 and rapamycin are known to inhibit FKBP11's peptidyl-prolyl isomerase activity.

### REFERENCES

- Galat, A. 1993. Peptidylproline *cis-trans*-isomerases: immunophilins. *Eur. J. Biochem.* 216: 689-707.
- Göthel, S.F., et al. 1999. Peptidyl-prolyl *cis-trans* isomerases, a superfamily of ubiquitous folding catalysts. *Cell. Mol. Life Sci.* 55: 423-436.
- Patterson, C.E., et al. 2002. Genomic organization of mouse and human 65 kDa FK506-binding protein genes and evolution of the FKBP multigene family. *Genomics* 79: 881-889.
- Breiman, A., et al. 2002. The involvement of mammalian and plant FK506-binding proteins (FKBPs) in development. *Transgenic Res.* 11: 321-335.
- Galat, A. 2003. Peptidylprolyl *cis-trans* isomerases (immunophilins): biological diversity—targets—functions. *Curr. Top. Med. Chem.* 3: 1315-1347.
- Rulten, S.L., et al. 2006. The human FK506-binding proteins: characterization of human FKBP19. *Mamm. Genome* 17: 322-331.
- Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2006. Johns Hopkins University, Baltimore, MD. MIM Number: 610571. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Kang, C.B., et al. 2008. FKBP family proteins: immunophilins with versatile biological functions. *Neurosignals* 16: 318-325.

### CHROMOSOMAL LOCATION

Genetic locus: FKBP11 (human) mapping to 12q13.12.

### PRODUCT

FKBP11 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see FKBP11 shRNA Plasmid (h): sc-75021-SH and FKBP11 shRNA (h) Lentiviral Particles: sc-75021-V as alternate gene silencing products.

For independent verification of FKBP11 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-75021A, sc-75021B and sc-75021C.

### STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

### APPLICATIONS

FKBP11 siRNA (h) is recommended for the inhibition of FKBP11 expression in human cells.

### SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

### GENE EXPRESSION MONITORING

FKBP11 (D-3): sc-398700 is recommended as a control antibody for monitoring of FKBP11 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker<sup>™</sup> Molecular Weight Standards: sc-2035, UltraCruz<sup>®</sup> Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz<sup>®</sup> Mounting Medium: sc-24941 or UltraCruz<sup>®</sup> Hard-set Mounting Medium: sc-359850.

### RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor FKBP11 gene expression knockdown using RT-PCR Primer: FKBP11 (h)-PR: sc-75021-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

### RESEARCH USE

For research use only, not for use in diagnostic procedures.

### PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.