

Gle1 siRNA (m): sc-75139

BACKGROUND

Protein transport across the nucleus is a selective, multi-step process involving several cytoplasmic factors that mediate protein passage through the nuclear pore complex (NPC). Gle1, also known as GLE1L, is a 698 amino acid protein that localizes to both the nucleus and the cytoplasm and belongs to the Gle1 family. Expressed as two alternatively spliced isoforms, Gle1 associates with the NPC and is required for the transport of poly(A)-containing mRNAs from the nucleus to the cytoplasm. Defects in the gene encoding Gle1 are the cause of lethal congenital contracture syndrome type 1 (LCCS1) and lethal arthrogryposis with anterior horn cell disease (LAAHD), the former of which is characterized by early fetal hydrops and akinesia, micrognathia, pulmonary hypoplasia, pterygia and prenatal death, while the latter is associated with respiratory failure.

REFERENCES

1. Watkins, J.L., et al. 1998. The human homologue of *Saccharomyces cerevisiae* Gle1p is required for poly(A)⁺ RNA export. *Proc. Natl. Acad. Sci. USA* 95: 6779-6784.
2. Kendirgi, F., et al. 2003. An essential role for hGle1 nucleocytoplasmic shuttling in mRNA export. *J. Cell Biol.* 160: 1029-1040.
3. Rayala, H.J., et al. 2004. The mRNA export factor human Gle1 interacts with the nuclear pore complex protein Nup155. *Mol. Cell Proteomics* 3: 145-155.
4. Kendirgi, F., et al. 2005. Interaction between the shuttling mRNA export factor Gle1 and the nucleoporin hCG1: a conserved mechanism in the export of HSP 70 mRNA. *Mol. Biol. Cell* 16: 4304-4315.
5. Kutay, U. and Panse, V.G. 2008. Gle1 does double duty. *Cell* 134: 564-566.
6. Bolger, T.A., et al. 2008. The mRNA export factor Gle1 and inositol hexakisphosphate regulate distinct stages of translation. *Cell* 134: 624-633.
7. Nousiainen, et al. 2008. Mutations in mRNA export mediator Gle1 result in a fetal motoneuron disease. *Nat. Genet.* 40: 155-157.
8. Online Mendelian Inheritance in Man, OMIM[™]. 2008. Johns Hopkins University, Baltimore, MD. MIM Number: 603371. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: Gle1 (mouse) mapping to 2 B.

PRODUCT

Gle1 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Gle1 shRNA Plasmid (m): sc-75139-SH and Gle1 shRNA (m) Lentiviral Particles: sc-75139-V as alternate gene silencing products.

For independent verification of Gle1 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-75139A, sc-75139B and sc-75139C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCL, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Gle1 siRNA (m) is recommended for the inhibition of Gle1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

Gle1 (G-9): sc-514796 is recommended as a control antibody for monitoring of Gle1 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Gle1 gene expression knockdown using RT-PCR Primer: Gle1 (m)-PR: sc-75139-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.