

# HORMAD1 siRNA (m): sc-75274

## BACKGROUND

HORMAD1 (HORMA domain containing 1), also known as NOHMA (newborn ovary HORMA protein) or CT46 (cancer/testis antigen 46), is a 394 amino acid testis-specific protein. HORMAD1 contains one HORMA (Hop1p, Rev7p and MAD2) domain (a domain involved in chromatin binding) that makes up the entire full-length sequence of the protein. Proteins with HORMA domains are typically involved in modulating chromatin dynamics and structure. The HORMA domain is believed to act as an adaptor, recruiting other proteins to chromatin states that result from nonattachment to the mitotic spindle or from DNA double-strand breaks. HORMAD1 is a putative meiotic protein, as is suggested by its 25.8% homology with the yeast protein Hop1 (a meiosis-specific protein). In addition, HORMAD1 is overexpressed in a variety of carcinomas, including breast, lung, endometrial, colon, bladder and esophageal cancers.

## REFERENCES

1. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 609824. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
2. Pangas, S.A., Yan, W., Matzuk, M.M. and Rajkovic, A. 2004. Restricted germ cell expression of a gene encoding a novel mammalian HORMA domain-containing protein. *Gene Expr. Patterns* 5: 257-263.
3. Chen, Y.T., Venditti, C.A., Theiler, G., Stevenson, B.J., Iseli, C., Gure, A.O., Jongeneel, C.V., Old, L.J. and Simpson, A.J. 2005. Identification of CT46/HORMAD1, an immunogenic cancer/testis antigen encoding a putative meiosis-related protein. *Cancer Immun.* 5: 9.
4. Aung, P.P., Que, N., Mitani, Y., Nakayama, H., Yoshida, K., Noguchi, T., Bosserhoff, A.K. and Yasui, W. 2006. Systematic search for gastric cancer-specific genes based on SAGE data: melanoma inhibitory activity and matrix metalloproteinase-10 are novel prognostic factors in patients with gastric cancer. *Oncogene* 25: 2546-2557.
5. Adelaide, J., Finetti, P., Bekhouche, I., Repellini, L., Geneix, J., Sircoulomb, F., Charafe-Jauffret, E., Cervera, N., Desplans, J., Parzy, D., Schoenmakers, E., Viens, P., Jacquemier, J., Birnbaum, D., Bertucci, F. and Chaffanet, M. 2007. Integrated profiling of basal and luminal breast cancers. *Cancer Res.* 67: 11565-11575.

## CHROMOSOMAL LOCATION

Genetic locus: Hormad1 (mouse) mapping to 3 F2.1.

## PRODUCT

HORMAD1 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see HORMAD1 shRNA Plasmid (m): sc-75274-SH and HORMAD1 shRNA (m) Lentiviral Particles: sc-75274-V as alternate gene silencing products.

For independent verification of HORMAD1 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-75274A, sc-75274B and sc-75274C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

HORMAD1 siRNA (m) is recommended for the inhibition of HORMAD1 expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor HORMAD1 gene expression knockdown using RT-PCR Primer: HORMAD1 (m)-PR: sc-75274-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.