

LSm14B siRNA (h): sc-75704

BACKGROUND

Sm and Sm-like (LSm) proteins form donut-shaped, ubiquitously expressed heptameric complexes that are involved in various steps of RNA metabolism, including RNA-protein interactions and structural changes that are required during ribosomal subunit assembly. LSm14B, also known as C20orf40, FAM61B or LSM13, is a 385 amino acid protein that exists as multiple alternatively spliced isoforms and may play a role in RNA-related events. The gene encoding LSm14B maps to human chromosome 20. Comprising approximately 2% of the human genome, chromosome 20 contains nearly 63 million bases that encode over 600 genes, some of which are associated with Creutzfeldt-Jakob disease, amyotrophic lateral sclerosis, spinal muscular atrophy, ring chromosome 20 epilepsy syndrome and Alagille syndrome.

REFERENCES

1. Wilusz, C.J. and Wilusz, J. 2005. Eukaryotic LSm proteins: lessons from bacteria. *Nat. Struct. Mol. Biol.* 12: 1031-1036.
2. Wolin, S.L. and Wurtmann, E.J. 2006. Molecular chaperones and quality control in non-coding RNA biogenesis. *Cold Spring Harb. Symp. Quant. Biol.* 71: 505-511.
3. Barbee, S.A. and Evans, T.C. 2006. The Sm proteins regulate germ cell specification during early *C. elegans* embryogenesis. *Dev. Biol.* 291: 132-143.
4. Ville, D., et al. 2006. Early pattern of epilepsy in the ring chromosome 20 syndrome. *Epilepsia* 47: 543-549.
5. Arluison, V., et al. 2006. Three-dimensional structures of fibrillar Sm proteins: Hfq and other Sm-like proteins. *J. Mol. Biol.* 356: 86-96.
6. Tharun, S. 2008. Chapter 4 roles of eukaryotic LSm proteins in the regulation of mRNA function. *Int. Rev. Cell Mol. Biol.* 272: 149-189.
7. Zaric, B.L. and Kambach, C. 2008. Reconstitution of recombinant human LSm complexes for biochemical, biophysical, and cell biological studies. *Meth. Enzymol.* 448: 57-74.

CHROMOSOMAL LOCATION

Genetic locus: LSM14B (human) mapping to 20q13.33.

PRODUCT

LSm14B siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see LSm14B shRNA Plasmid (h): sc-75704-SH and LSm14B shRNA (h) Lentiviral Particles: sc-75704-V as alternate gene silencing products.

For independent verification of LSm14B (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-75704A, sc-75704B and sc-75704C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCL, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

LSm14B siRNA (h) is recommended for the inhibition of LSm14B expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor LSm14B gene expression knockdown using RT-PCR Primer: LSm14B (h)-PR: sc-75704-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.