

Maltase-glucoamylase siRNA (h): sc-75740

BACKGROUND

Maltase-glucoamylase, also known as MGAM, MG or MGA, is a 1,857 amino acid multi-pass membrane protein that localizes to the apical cell membrane and contains two P-type domains. Expressed in kidney, small intestine and granulocytes, Maltase-glucoamylase exists as a monomer that is thought to participate in an alternate pathway of starch digestion, specifically when luminal α -amylase activity is reduced because of immaturity or malnutrition. Maltase-glucoamylase is subject to posttranslational N- and O-glycosylation, as well as sulfation. The gene encoding Maltase-glucoamylase maps to human chromosome 7q34, which houses over 1,000 genes and comprises nearly 5% of the human genome. Defects in some of the genes localized to chromosome 7 have been linked to osteogenesis imperfecta, Williams-Beuren syndrome, Pendred syndrome, lissencephaly, citrullinemia and Shwachman-Diamond syndrome.

REFERENCES

1. Danielsen, E.M. 1987. Tyrosine sulfation, a post-translational modification of microvillar enzymes in the small intestinal enterocyte. *EMBO J.* 6: 2891-2896.
2. Naim, H.Y., Sterchi, E.E. and Lentze, M.J. 1988. Structure, biosynthesis, and glycosylation of human small intestinal Maltase-glucoamylase. *J. Biol. Chem.* 263: 19709-19717.
3. Nichols, B.L., Eldering, J., Avery, S., Hahn, D., Quaroni, A. and Sterchi, E. 1998. Human small intestinal Maltase-glucoamylase cDNA cloning. Homology to sucrase-isomaltase. *J. Biol. Chem.* 273: 3076-3081.
4. Nichols, B.L., Avery, S., Sen, P., Swallow, D.M., Hahn, D. and Sterchi, E. 2003. The Maltase-glucoamylase gene: common ancestry to sucrase-isomaltase with complementary starch digestion activities. *Proc. Natl. Acad. Sci. USA* 100: 1432-1437.
5. Online Mendelian Inheritance in Man, OMIM™. 2003. Johns Hopkins University, Baltimore, MD. MIM Number: 154360. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Ao, Z., Quezada-Calvillo, R., Sim, L., Nichols, B.L., Rose, D.R., Sterchi, E.E. and Hamaker, B.R. 2007. Evidence of native starch degradation with human small intestinal Maltase-glucoamylase (recombinant). *FEBS Lett.* 581: 2381-2388.
7. Tu, S., Kutlu, T., Pekin, O., Baglam, E., Kiyak, H. and Oral, O. 2008. Clinical, endocrine, and metabolic effects of acarbose, an α -glucosidase inhibitor, in overweight and nonoverweight patients with polycystic ovarian syndrome. *Fertil. Steril.* 90: 1144-1148.
8. Sim, L., Quezada-Calvillo, R., Sterchi, E.E., Nichols, B.L. and Rose, D.R. 2008. Human intestinal Maltase-glucoamylase: crystal structure of the N-terminal catalytic subunit and basis of inhibition and substrate specificity. *J. Mol. Biol.* 375: 782-792.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: MGAM (human) mapping to 7q34.

PRODUCT

Maltase-glucoamylase siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Maltase-glucoamylase shRNA Plasmid (h): sc-75740-SH and Maltase-glucoamylase shRNA (h) Lentiviral Particles: sc-75740-V as alternate gene silencing products.

For independent verification of Maltase-glucoamylase (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-75740A, sc-75740B and sc-75740C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Maltase-glucoamylase siRNA (h) is recommended for the inhibition of Maltase-glucoamylase expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Maltase-glucoamylase gene expression knockdown using RT-PCR Primer: Maltase-glucoamylase (h)-PR: sc-75740-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.