

Rab 18 siRNA (h): sc-76316

BACKGROUND

The Ras-related superfamily of guanine nucleotide binding proteins includes the R-Ras, Rap, Ral/Rec and Rho/Rab subfamilies all of which are thought to play an important role in either endocytosis or in biosynthetic protein transport. The process of transporting newly synthesized proteins from the endoplasmic reticulum (ER) to various stacks of the Golgi complex and to secretory vesicles involves the movement of carrier vesicles and requires Rab protein function. Rab proteins are also an integral part of endocytic pathways. Rab 18 is a 206 amino acid protein that is lipid-anchored to the cytoplasmic side of the cell membrane and belongs to the Rab subfamily of small GTPases. Expressed ubiquitously, Rab 18 plays a role in apical recycling and endocytosis and is also thought to be involved in protein transport between early endosomes and the plasma membrane.

REFERENCES

- McMurtrie, E.B., et al. 1997. Rab 17 and Rab 18, small GTPases with specificity for polarized epithelial cells: genetic mapping in the mouse. *Genomics* 45: 623-625.
- Schäfer, U., et al. 2000. Isolation and characterisation of the human Rab 18 gene after stimulation of endothelial cells with histamine. *FEBS Lett.* 466: 148-154.
- Online Mendelian Inheritance in Man, OMIM[™]. 2001. Johns Hopkins University, Baltimore, MD. MIM Number: 602207. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Dou, T., et al. 2005. Cloning and characterization of a novel splice variant of human Rab 18 gene (RAB18). *DNA Seq.* 16: 230-234.
- Martin, S., et al. 2005. Regulated localization of Rab 18 to lipid droplets: effects of lipolytic stimulation and inhibition of lipid droplet catabolism. *J. Biol. Chem.* 280: 42325-42335.
- Ozeki, S., et al. 2005. Rab 18 localizes to lipid droplets and induces their close apposition to the endoplasmic reticulum-derived membrane. *J. Cell Sci.* 118: 2601-2611.
- Dejgaard, S.Y., et al. 2008. Rab 18 and Rab 43 have key roles in ER-Golgi trafficking. *J. Cell Sci.* 121: 2768-2781.

CHROMOSOMAL LOCATION

Genetic locus: RAB18 (human) mapping to 10p12.1.

PRODUCT

Rab 18 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Rab 18 shRNA Plasmid (h): sc-76316-SH and Rab 18 shRNA (h) Lentiviral Particles: sc-76316-V as alternate gene silencing products.

For independent verification of Rab 18 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-76316A, sc-76316B and sc-76316C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Rab 18 siRNA (h) is recommended for the inhibition of Rab 18 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

Rab 18 (D-5): sc-393168 is recommended as a control antibody for monitoring of Rab 18 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Rab 18 gene expression knockdown using RT-PCR Primer: Rab 18 (h)-PR: sc-76316-PR (20 μ l, 415 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.