

# SerpinB8 siRNA (h): sc-76478

## BACKGROUND

The serine proteinase inhibitors (serpins) comprise a superfamily of proteins with a diverse set of functions, including the control of blood coagulation, complement activation, programmed cell death and tissue development. SerpinB8 (serpin peptidase inhibitor, clade B (ovalbumin), member 8), also known as P18 or CAP2, is a 374 amino acid protein that localizes to the cytoplasm and belongs to the serine proteinase inhibitor family. Expressed at high levels in lung, liver, heart and skeletal muscle, SerpinB8, which exists as multiple alternatively spliced isoforms, is thought to play a role in platelet aggregation and platelet-regulated pathophysiological responses. The gene encoding SerpinB8 maps to a cluster of Serpin genes on human chromosome 18, which houses over 300 protein-coding genes and contains nearly 76 million bases.

## REFERENCES

- Huber, R. and Carrell, R.W. 1989. Implications of the three-dimensional structure of  $\alpha$  1-antitrypsin for structure and function of serpins. *Biochemistry* 28: 8951-8966.
- Sprecher, C.A., et al. 1995. Molecular cloning, expression, and partial characterization of two novel members of the ovalbumin family of serine proteinase inhibitors. *J. Biol. Chem.* 270: 29854-29861.
- Bartuski, A.J., et al. 1997. Cytoplasmic antiproteinase 2 (PI8) and bomapin (PI10) map to the serpin cluster at 18q21.3. *Genomics* 43: 321-328.
- Dahlen, J.R., et al. 1998. Inhibition of soluble recombinant furin by human proteinase inhibitor 8. *J. Biol. Chem.* 273: 1851-1854.
- Strik, M.C., et al. 2002. Distribution of the human intracellular serpin protease inhibitor 8 in human tissues. *J. Histochem. Cytochem.* 50: 1443-1454.
- Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 601697. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Pike, R.N., et al. 2005. Control of the coagulation system by serpins. Getting by with a little help from glycosaminoglycans. *FEBS J.* 272: 4842-4851.

## CHROMOSOMAL LOCATION

Genetic locus: SERPINB8 (human) mapping to 18q22.1.

## PRODUCT

SerpinB8 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see SerpinB8 shRNA Plasmid (h): sc-76478-SH and SerpinB8 shRNA (h) Lentiviral Particles: sc-76478-V as alternate gene silencing products.

For independent verification of SerpinB8 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-76478A, sc-76478B and sc-76478C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

SerpinB8 siRNA (h) is recommended for the inhibition of SerpinB8 expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

SerpinB8 (PI-8): sc-101371 is recommended as a control antibody for monitoring of SerpinB8 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker<sup>™</sup> Molecular Weight Standards: sc-2035, UltraCruz<sup>®</sup> Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz<sup>®</sup> Mounting Medium: sc-24941 or UltraCruz<sup>®</sup> Hard-set Mounting Medium: sc-359850.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor SerpinB8 gene expression knockdown using RT-PCR Primer: SerpinB8 (h)-PR: sc-76478-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.