

# cyclin O siRNA (m): sc-76803

## BACKGROUND

Cell proliferation is controlled at specific stages of the cell cycle by distinct protein kinase complexes. These complexes consist of a catalytic subunit associating with a specific regulatory subunit to form the active kinase. The cyclins, which include cyclin A, B, C, D, E, F, G, H, I, K, L, O, T and their related proteins, comprise the regulatory subunits of these kinase complexes. The controlled activation of the kinase complexes at various intervals of the cell cycle is regulated by the availability of the cyclins to the catalytic subunit. Unlike the catalytic subunit, which is expressed continually, the expression and stability of the regulatory subunit fluctuates depending on the stage of the cell cycle, thereby regulating kinase activity. Cyclin O, also known as CCNO, is a 350 amino acid protein that belongs to the cyclin family and is encoded by a gene located on human chromosome 5q11.2. Cyclin O may play an important role in oocyte meiotic cell cycle.

## REFERENCES

- Gallant, P. and Nigg, E.A. 1994. Identification of a novel vertebrate cyclin: cyclin B3 shares properties with both A- and B-type cyclins. *EMBO J.* 13: 595-605.
- Mikulits, W., Dolznig, H., Edelmann, H., Sauer, T., Deiner, E.M., Ballou, L., Beug, H. and Müllner, E.W. 1997. Dynamics of cell cycle regulators: artifact-free analysis by recultivation of cells synchronized by centrifugal elutriation. *DNA Cell Biol.* 16: 849-859.
- Kolonin, M.G. and Finley, R.L. 2000. A role for cyclin J in the rapid nuclear division cycles of early *Drosophila* embryogenesis. *Dev. Biol.* 227: 661-672.
- Kong, M., Barnes, E.A., Ollendorff, V. and Donoghue, D.J. 2000. Cyclin F regulates the nuclear localization of cyclin B1 through a cyclin-cyclin interaction. *EMBO J.* 19: 1378-1388.
- Wikman, H. and Kettunen, E. 2006. Regulation of the G<sub>1</sub>/S phase of the cell cycle and alterations in the RB pathway in human lung cancer. *Expert Rev. Anticancer Ther.* 6: 515-530.
- Roig, M.B., Roset, R., Ortet, L., Balsiger, N.A., Anfosso, A., Cabellos, L., Garrido, M., Alameda, F., Brady, H.J. and Gil-Gómez, G. 2009. Identification of a novel cyclin required for the intrinsic apoptosis pathway in lymphoid cells. *Cell Death Differ.* 16: 230-243.
- Roset, R. and Gil-Gómez, G. 2009. Measurement of changes in Cdk2 and cyclin O-associated kinase activity in apoptosis. *Methods Mol. Biol.* 559: 161-172.
- Ma, J.Y., Ou-Yang, Y.C., Luo, Y.B., Wang, Z.B., Hou, Y., Han, Z.M., Liu, Z., Schatten, H. and Sun, Q.Y. 2013. Cyclin O regulates germinal vesicle breakdown in mouse oocytes. *Biol. Reprod.* 88: 110.

## CHROMOSOMAL LOCATION

Genetic locus: Ccno (mouse) mapping to 13 D2.2.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

## PRODUCT

cyclin O siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see cyclin O shRNA Plasmid (m): sc-76803-SH and cyclin O shRNA (m) Lentiviral Particles: sc-76803-V as alternate gene silencing products.

For independent verification of cyclin O (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-76803A, sc-76803B and sc-76803C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

cyclin O siRNA (m) is recommended for the inhibition of cyclin O expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor cyclin O gene expression knockdown using RT-PCR Primer: cyclin O (m)-PR: sc-76803-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.