



ZFHX4 siRNA (h): sc-76956

BACKGROUND

Zinc-finger proteins contain DNA-binding domains and have a wide variety of functions, most of which encompass some form of transcriptional activation or repression. The majority of zinc-finger proteins contain a Krüppel-type DNA binding domain and a KRAB domain, which are thought to interact with KAP1, thereby recruiting histone modifying proteins. ZFHX4, also designated ZFH4 or ZHF4, is a 3,567 amino acid nuclear protein expressed in brain, skeletal muscle and liver with very low expression in stomach. Belonging to the Krüppel C₂H₂-type zinc-finger protein family, ZFHX4 may play a role in neural and muscle differentiation and may be involved in transcriptional regulation. ZFHX4 contains 20 C₂H₂-type zinc fingers and 4 homeobox DNA-binding domains, and exists as 3 alternatively spliced variants.

REFERENCES

1. DeLisi, C., et al. 1981. Effect of nonspecific forces and finite receptor number on rate constants of ligand—cell bound-receptor interactions. *Proc. Natl. Acad. Sci. USA* 78: 5569-5572.
2. Skapek, S.X., et al. 2000. Cloning and characterization of a novel Krüppel-associated box family transcriptional repressor that interacts with the retinoblastoma gene product, RB. *J. Biol. Chem.* 275: 7212-7223.
3. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 606940. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
4. Nakamura, M., et al. 2004. A novel subfamily of zinc finger genes involved in embryonic development. *J. Cell. Biochem.* 93: 887-895.
5. Nakashima, M., et al. 2008. Genome-wide linkage analysis and mutation analysis of hereditary congenital blepharoptosis in a Japanese family. *J. Hum. Genet.* 53: 34-41.
6. Comabella, M., et al. 2009. Genome-wide scan of 500,000 single-nucleotide polymorphisms among responders and nonresponders to interferon β therapy in multiple sclerosis. *Arch. Neurol.* 66: 972-978.
7. Thomas, J.H., et al. 2009. Evolution of C₂H₂-zinc finger genes revisited. *BMC Evol. Biol.* 9: 51.

CHROMOSOMAL LOCATION

Genetic locus: ZFHX4 (human) mapping to 8q21.11.

PRODUCT

ZFHX4 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ZFHX4 shRNA Plasmid (h): sc-76956-SH and ZFHX4 shRNA (h) Lentiviral Particles: sc-76956-V as alternate gene silencing products.

For independent verification of ZFHX4 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-76956A, sc-76956B and sc-76956C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

ZFHX4 siRNA (h) is recommended for the inhibition of ZFHX4 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ZFHX4 gene expression knockdown using RT-PCR Primer: ZFHX4 (h)-PR: sc-76956-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.