

# cyclin M2 siRNA (m): sc-77062

## BACKGROUND

Cyclin M2, also known as CNM2 or ACDP2 (ancient conserved domain-containing protein 2), is an 875 amino acid multi-pass membrane protein that contains two CBS domains and belongs to the ACDP family. Expressed in a variety of tissues with highest expression in placenta, brain and kidney, cyclin M2 functions as a divalent metal cation transporter that mediates the transport of several different metal cations, including  $Mg^{2+}$ ,  $Co^{2+}$  and  $Fe^{2+}$ . Cyclin M2 exists as multiple alternatively spliced isoforms and, contrary to its name, exhibits no cyclin-like function *in vivo*. The gene encoding cyclin M2 maps to human chromosome 10, which houses over 1,200 genes and comprises nearly 4.5% of the human genome. Defects in some of the genes that map to chromosome 10 are associated with Charcot-Marie-Tooth disease, Jackson-Weiss syndrome, Usher syndrome, nonsyndromic deafness, Wolman's syndrome, Cowden syndrome, multiple endocrine neoplasia type 2 and porphyria.

## REFERENCES

- Berger, P., Young, P. and Suter, U. 2002. Molecular cell biology of Charcot-Marie-Tooth disease. *Neurogenetics* 4: 1-15.
- Wang, C.Y., Shi, J.D., Yang, P., Kumar, P.G., Li, Q.Z., Run, Q.G., Su, Y.C., Scott, H.S., Kao, K.J. and She, J.X. 2003. Molecular cloning and characterization of a novel gene family of four ancient conserved domain proteins (ACDP). *Gene* 306: 37-44.
- Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2003. Johns Hopkins University, Baltimore, MD. MIM Number: 607803. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
- Wang, C.Y., Yang, P., Shi, J.D., Purohit, S., Guo, D., An, H., Gu, J.G., Ling, J., Dong, Z. and She, J.X. 2004. Molecular cloning and characterization of the mouse Acdp gene family. *BMC Genomics* 5: 7.
- Goytain, A. and Quamme, G.A. 2005. Functional characterization of ACDP2 (ancient conserved domain protein), a divalent metal transporter. *Physiol. Genomics* 22: 382-389.
- Schmitz, C., Deason, F. and Perraud, A.L. 2007. Molecular components of vertebrate  $Mg^{2+}$ -homeostasis regulation. *Magnesium Res.* 20: 6-18.
- Sontia, B. and Touyz, R.M. 2007. Magnesium transport in hypertension. *Pathophysiology* 14: 205-211.

## CHROMOSOMAL LOCATION

Genetic locus: Cnm2 (mouse) mapping to 19 C3.

## PRODUCT

cyclin M2 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see cyclin M2 shRNA Plasmid (m): sc-77062-SH and cyclin M2 shRNA (m) Lentiviral Particles: sc-77062-V as alternate gene silencing products.

For independent verification of cyclin M2 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-77062A, sc-77062B and sc-77062C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

cyclin M2 siRNA (m) is recommended for the inhibition of cyclin M2 expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor cyclin M2 gene expression knockdown using RT-PCR Primer: cyclin M2 (m)-PR: sc-77062-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.