

cystatin 11 siRNA (h): sc-77080

BACKGROUND

The cystatins comprise a superfamily of proteins that contain multiple cystatin-like sequences and may function as active cysteine protease inhibitors. Cystatin proteins are divided into three distinct subfamilies, designated type 1 cystatins (stefins), type 2 cystatins and kininogens. The type 2 cystatins are expressed in human fluids and secretions and are encoded by genes that localize to a cystatin locus on chromosome 20. Cystatin 11, also known as SC13 or CST8L, is a 138 amino acid protein that localizes to both the nucleus and the cytoplasm and belongs to the type 2 subfamily of cystatins. Detected in epididymal and epithelial tissue, cystatin 11 functions as an antibacterial protein that is thought to play a role in sperm maturation and fertilization. Multiple isoforms of cystatin 11 exist due to alternative splicing events.

REFERENCES

1. Saitoh, E., Isemura, S., Sanada, K., Kim, H.S., Smithies, O. and Maeda, N. 1988. Cystatin superfamily. Evidence that family II cystatin genes are evolutionarily related to family III cystatin genes. *Biol. Chem. Hoppe Seyler* 369: 191-197.
2. Dickinson, D.P., Thiesse, M., Dempsey, L.D. and Millar, S.J. 1993. Genomic cloning, physical mapping, and expression of human type 2 cystatin genes. *Crit. Rev. Oral Biol. Med.* 4: 573-580.
3. Thiesse, M., Millar, S.J. and Dickinson, D.P. 1994. The human type 2 cystatin gene family consists of eight to nine members, with at least seven genes clustered at a single locus on human chromosome 20. *DNA Cell Biol.* 13: 97-116.
4. Dickinson, D.P., Zhao, Y., Thiesse, M. and Siciliano, M.J. 1994. Direct mapping of seven genes encoding human type 2 cystatins to a single site located at 20p11.2. *Genomics* 24: 172-175.
5. Brown, W.M. and Dziegielewska, K.M. 1997. Friends and relations of the cystatin superfamily—new members and their evolution. *Protein Sci.* 6: 5-12.
6. Hamil, K.G., Liu, Q., Sivashanmugam, P., Yenugu, S., Soundararajan, R., Grossman, G., Richardson, R.T., Zhang, Y.L., O'Rand, M.G., Petrusz, P., French, F.S. and Hall, S.H. 2002. Cystatin 11: a new member of the cystatin type 2 family. *Endocrinology* 143: 2787-2796.
7. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 609731. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
8. Cornwall, G.A. and Hsia, N. 2003. A new subgroup of the family 2 cystatins. *Mol. Cell. Endocrinol.* 200: 1-8.
9. Frygeli, J., Oscarson, M., Nordqvist, K., Wedell, A. and Tökönen, V. 2007. The reproductive tissue specific cystatin subgroup of genes: expression during gonadal development in wildtype and testatin knockout animals. *Sex Dev.* 1: 363-372.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: CST11 (human) mapping to 20p11.21.

PRODUCT

cystatin 11 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see cystatin 11 shRNA Plasmid (h): sc-77080-SH and cystatin 11 shRNA (h) Lentiviral Particles: sc-77080-V as alternate gene silencing products.

For independent verification of cystatin 11 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-77080A, sc-77080B and sc-77080C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

cystatin 11 siRNA (h) is recommended for the inhibition of cystatin 11 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RESEARCH USE

For research use only, not for use in diagnostic procedures.