



DPPA5 siRNA (h): sc-77180

BACKGROUND

DPPA5 (developmental pluripotency associated 5), also known as ESG1 (embryonal stem cell-specific gene 1), is a 116 amino acid protein that localizes to the cytoplasm and contains one KH domain. Expressed in embryonic germ (EG), primordial germ (PG) and embryonic stem (ES) cells, DPPA5 plays an important role in the maintenance of ES cell pluripotency and may be necessary for proper embryogenesis. The gene encoding DPPA5 maps to human chromosome 6, which contains 170 million base pairs and comprises nearly 6% of the human genome. Deletion of a portion of the q arm of chromosome 6 is associated with early onset intestinal cancer, suggesting the presence of a cancer susceptibility locus. Additionally, Porphyria cutanea tarda, Parkinson's disease, Stickler syndrome and a susceptibility to bipolar disorder are all associated with genes that map to chromosome 6.

REFERENCES

1. Kim, S.K., Suh, M.R., Yoon, H.S., Lee, J.B., Oh, S.K., Moon, S.Y., Moon, S.H., Lee, J.Y., Hwang, J.H., Cho, W.J. and Kim, K.S. 2005. Identification of developmental pluripotency associated 5 expression in human pluripotent stem cells. *Stem Cells* 23: 458-462.
2. Amano, H., Itakura, K., Maruyama, M., Ichisaka, T., Nakagawa, M. and Yamanaka, S. 2006. Identification and targeted disruption of the mouse gene encoding ESG1 (PH34/ECAT2/DPPA5). *BMC Dev. Biol.* 6: 11.
3. Lagarkova, M.A., Volchkov, P.Y., Lyakisheva, A.V., Philonenko, E.S. and Kiselev, S.L. 2006. Diverse epigenetic profile of novel human embryonic stem cell lines. *Cell Cycle* 5: 416-420.
4. Pierre, A., Gautier, M., Callebaut, I., Bontoux, M., Jeanpierre, E., Pontarotti, P. and Monget, P. 2007. Atypical structure and phylogenomic evolution of the new eutherian oocyte- and embryo-expressed KHDC1/DPPA5/ECAT1/OOEP gene family. *Genomics* 90: 583-594.
5. Online Mendelian Inheritance in Man, OMIM™. 2007. Johns Hopkins University, Baltimore, MD. MIM Number: 611111. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Izadyar, F., Pau, F., Marh, J., Slepko, N., Wang, T., Gonzalez, R., Ramos, T., Howerton, K., Sayre, C. and Silva, F. 2008. Generation of multipotent cell lines from a distinct population of male germ line stem cells. *Reproduction* 135: 771-784.

CHROMOSOMAL LOCATION

Genetic locus: DPPA5 (human) mapping to 6q13.

PRODUCT

DPPA5 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see DPPA5 shRNA Plasmid (h): sc-77180-SH and DPPA5 shRNA (h) Lentiviral Particles: sc-77180-V as alternate gene silencing products.

For independent verification of DPPA5 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-77180A, sc-77180B and sc-77180C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

DPPA5 siRNA (h) is recommended for the inhibition of DPPA5 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor DPPA5 gene expression knockdown using RT-PCR Primer: DPPA5 (h)-PR: sc-77180-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.