

EF-1 β siRNA (h): sc-77236

BACKGROUND

EF-1 (elongation factor-1) is a multi-protein complex that is responsible for the delivery of aminoacyl-tRNAs to the ribosome. The EF-1 protein complex is composed of four subunits: EF-1 α , EF-1 β , EF-1 δ and EF-1 γ . EF-1 β is a 225 amino acid member of the EF-1 β /EF-1 δ protein family and contains one GST C-terminal domain. Both EF-1 β and EF-1 δ stimulate the exchange of GDP bound to EF-1 α to GTP. Phosphorylation of EF-1 β affects the GDP/GTP exchange rate. Alternative splicing events of the gene that encodes EF-1 β result in three transcript variants that differ only in the 5' UTR.

REFERENCES

1. Cho, D.I., et al. 2003. Direct and biochemical interaction between dopamine D3 receptor and elongation factor-1 β γ . *Life Sci.* 73: 2991-3004.
2. Cans, C., et al. 2003. Translationally controlled tumor protein acts as a guanine nucleotide dissociation inhibitor on the translation elongation factor eEF1A. *Proc. Natl. Acad. Sci. USA* 100: 13892-13897.
3. Ito, T., et al. 2004. Solution structure of human initiation factor eIF2 α reveals homology to the elongation factor eEF1B. *Structure* 12: 1693-1704.
4. McCracken, S., et al. 2005. Proteomic analysis of SRm160-containing complexes reveals a conserved association with cohesin. *J. Biol. Chem.* 280: 42227-42236.
5. Le Sourd, F., et al. 2006. eEF1B: at the dawn of the 21st century. *Biochim. Biophys. Acta* 1759: 13-31.
6. Le Sourd, F., et al. 2006. Cellular coexistence of two high molecular subsets of eEF1B complex. *FEBS Lett.* 580: 2755-2760.
7. Deng, S.S., et al. 2006. Comparative proteome analysis of breast cancer and adjacent normal breast tissues in human. *Genomics Proteomics Bioinformatics* 4: 165-172.
8. Mazan-Mamczarz, K., et al. 2006. Translational repression by RNA-binding protein TIAR. *Mol. Cell. Biol.* 26: 2716-2727.
9. Byun, H.O., et al. 2009. Cathepsin D and eukaryotic translation elongation factor 1 as promising markers of cellular senescence. *Cancer Res.* 69: 4638-4647.

CHROMOSOMAL LOCATION

Genetic locus: EEF1B2 (human) mapping to 2q33.3.

PRODUCT

EF-1 β siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see EF-1 β shRNA Plasmid (h): sc-77236-SH and EF-1 β shRNA (h) Lentiviral Particles: sc-77236-V as alternate gene silencing products.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

EF-1 β siRNA (h) is recommended for the inhibition of EF-1 β expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

EF-1 β (3A5): sc-517178 is recommended as a control antibody for monitoring of EF-1 β gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor EF-1 β gene expression knockdown using RT-PCR Primer: EF-1 β (h)-PR: sc-77236-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.