

GCP3 siRNA (m): sc-77349

BACKGROUND

The γ -Tubulin complex is composed of γ Tubulin and the γ -Tubulin complex-associated proteins GCP2, GCP3, GCP4, GCP5 and GCP6, all of which are essential components of microtubule organizing centers. γ -Tubulin complex components are localized to both the centrosome, where they are involved in microtubule nucleation, and to the cytoplasm, where they exist as soluble complexes that can be recruited to the centrosome as needed. Although the GCP proteins are related, they have distinct roles which contribute to the proper function of the γ -Tubulin complex. GCP3 (γ -Tubulin complex component 3), also known as TUBGCP3 or SPBC98, localizes to the centrosome and is a ubiquitously expressed 907 amino acid member of the γ -Tubulin complex. Like GCP2 and γ Tubulin, GCP3 is conserved in all eukaryotes, suggesting that it is part of a core unit involved in eukaryotic microtubule nucleation. Three isoforms of GCP3 exist due to alternative splicing events.

REFERENCES

1. Murphy, S.M., Urbani, L. and Stearns, T. 1998. The mammalian γ Tubulin complex contains homologues of the yeast spindle pole body components spc97p and spc98p. *J. Cell Biol.* 141: 663-674.
2. Fava, F., Raynaud-Messina, B., Leung-Tack, J., Mazzolini, L., Li, M., Guillemot, J.C., Cachot, D., Tollon, Y., Ferrara, P. and Wright, M. 1999. Human 76p: a new member of the γ Tubulin-associated protein family. *J. Cell Biol.* 147: 857-868.
3. Hillman, R.T., Green, R.E. and Brenner, S.E. 2004. An unappreciated role for RNA surveillance. *Genome Biol.* 5: R8.
4. Rush, J., Moritz, A., Lee, K.A., Guo, A., Goss, V.L., Spek, E.J., Zhang, H., Zha, X.M., Polakiewicz, R.D. and Comb, M.J. 2005. Immunoaffinity profiling of tyrosine phosphorylation in cancer cells. *Nat. Biotechnol.* 23: 94-101.
5. Delgehyr, N., Sillibourne, J. and Bornens, M. 2005. Microtubule nucleation and anchoring at the centrosome are independent processes linked by Ninein function. *J. Cell Sci.* 118: 1565-1575.
6. Arbildua, J.J., Brunet, J.E., Jameson, D.M., López, M., Nova, E., Lagos, R. and Monasterio, O. 2006. Fluorescence resonance energy transfer and molecular modeling studies on 4',6-diamidino-2-phenylindole (DAPI) complexes with Tubulin. *Protein Sci.* 15: 410-419.
7. Stirling, P.C., Cuellar, J., Alfaro, G.A., El Khadali, F., Beh, C.T., Valpuesta, J.M., Melki, R. and Leroux, M.R. 2006. PhLP3 modulates CCT-mediated Actin and Tubulin folding via ternary complexes with substrates. *J. Biol. Chem.* 281: 7012-7021.
8. Van Thuan, N., Wakayama, S., Kishigami, S. and Wakayama, T. 2006. Donor centrosome regulation of initial spindle formation in mouse somatic cell nuclear transfer: roles of γ Tubulin and nuclear mitotic apparatus protein 1. *Biol. Reprod.* 74: 777-787.
9. Lüders, J., Patel, U.K. and Stearns, T. 2006. GCP-WD is a γ Tubulin targeting factor required for centrosomal and chromatin-mediated microtubule nucleation. *Nat. Cell Biol.* 8: 137-147.

CHROMOSOMAL LOCATION

Genetic locus: *Tubgcp3* (mouse) mapping to 8 A1.1.

PRODUCT

GCP3 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see GCP3 shRNA Plasmid (m): sc-77349-SH and GCP3 shRNA (m) Lentiviral Particles: sc-77349-V as alternate gene silencing products.

For independent verification of GCP3 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-77349A, sc-77349B and sc-77349C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

GCP3 siRNA (m) is recommended for the inhibition of GCP3 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor GCP3 gene expression knockdown using RT-PCR Primer: GCP3 (m)-PR: sc-77349-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.