

DDX27 siRNA (m): sc-77412

BACKGROUND

DDX27 (probable ATP-dependent RNA helicase DDX27) is a 796 amino acid protein encoded by the human gene DDX27. This protein belongs to the DEAD box helicase family, DDX27/DRS1 subfamily and contains one helicase ATP-binding domain and one helicase C-terminal domain. DDX27 is a nuclear protein found on chromosome 5 that likely functions as an ATP-dependent RNA helicase. RNA helicases are highly conserved enzymes that utilize the energy derived from NTP hydrolysis to modulate the structure of RNA. RNA Helicases participate in all biological processes that involve RNA, including transcription, splicing and translation.

REFERENCES

1. O'Day, C.L., Chavanikamannil, F. and Abelson, J. 1996. 18S rRNA processing requires the RNA helicase-like protein Rrp3. *Nucleic Acids Res.* 24: 3201-3207.
2. Doorbar, J., Elston, R.C., Napthine, S., Raj, K., Medcalf, E., Jackson, D., Coleman, N., Griffin, H.M., Masterson, P., Stacey, S., Mengistu, Y. and Dunlop, J. 2000. The E1E4 protein of human papillomavirus type 16 associates with a putative RNA helicase through sequences in its C terminus. *J. Virol.* 74: 10081-10095.
3. Regard, J.B., Scheek, S., Borbiev, T., Lanahan, A.A., Schneider, A., Demetriades, A.M., Hiemisch, H., Barnes, C.A., Verin, A.D. and Worley, P.F. 2004. Verge: a novel vascular early response gene. *J. Neurosci.* 24: 4092-4103.
4. Abdelhaleem, M. 2005. RNA Helicases: regulators of differentiation. *Clin. Biochem.* 38: 499-503.
5. Lee, J.H., Rho, S.B. and Chun, T. 2005. GABA_A receptor-associated protein (GABARAP) induces apoptosis by interacting with DEAD (Asp-Glu-Ala-Asp/His) box polypeptide 47 (DDX 47). *Biotechnol. Lett.* 27: 623-628.
6. Nousiainen, M., Sillje, H.H., Sauer, G., Nigg, E.A. and Körner, R. 2006. Phosphoproteome analysis of the human mitotic spindle. *Proc. Natl. Acad. Sci. USA* 103: 5391-5396.
7. Sekiguchi, T., Hayano, T., Yanagida, M., Takahashi, N. and Nishimoto, T. 2006. NOP132 is required for proper nucleolus localization of DEAD-box RNA helicase DDX47. *Nucleic Acids Res.* 34: 4593-4608.

CHROMOSOMAL LOCATION

Genetic locus: Ddx27 (mouse) mapping to 2 H3.

PRODUCT

DDX27 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see DDX27 shRNA Plasmid (m): sc-77412-SH and DDX27 shRNA (m) Lentiviral Particles: sc-77412-V as alternate gene silencing products.

For independent verification of DDX27 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-77412A, sc-77412B and sc-77412C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

DDX27 siRNA (m) is recommended for the inhibition of DDX27 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

DDX27 (2251C2a): sc-81074 is recommended as a control antibody for monitoring of DDX27 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor DDX27 gene expression knockdown using RT-PCR Primer: DDX27 (m)-PR: sc-77412-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.