



HR siRNA (h): sc-77654

BACKGROUND

HR (protein Hairless) is a 1,189 amino acid protein which is expressed as two isoforms produced by alternative splicing. The two isoforms are expressed in a variety of tissues in varying concentrations. Isoform 1 is more abundant than isoform 2 and is expressed at low levels in kidneys and testis, while isoform 2 is expressed abundantly in skin. Both isoforms are also present together in many tissues and are expressed strongly in small intestine and brain and weakly in trachea. HR is thought to be a transcription factor involved in hair growth. Hair growth occurs in three phases known as anagen, catagen and telogen, which are phases where growth, regression and rest, respectively, are taking place. By unknown mechanisms, HR is thought to regulate one of the hair growth phases and to work with vitamin D receptor (VDR) to regulate hair follicle cycling. Defects in HR may cause two serious ailments, known as alopecia universalis congenita (ALUNC) and atrichia with papular lesions (APL), which is also referred to as congenital atrichia. Both are autosomally recessive impairments. ALUNC is a rare condition in which hair follicles are produced without hair, while APL is a serious disease in which papillary lesions may cover the body and little to no hair is grown.

REFERENCES

1. Potter, G.B., et al. 2001. The Hairless gene mutated in congenital hair loss disorders encodes a novel nuclear receptor corepressor. *Genes Dev.* 15: 2687-2701.
2. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 602302. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
3. Djabali, K. and Christiano, A.M. 2004. Hairless contains a novel nuclear matrix targeting signal and associates with histone deacetylase 3 in nuclear speckles. *Differentiation* 72: 410-418.
4. Bergman, R., et al. 2005. The alopecias associated with vitamin D-dependent rickets type IIA and with Hairless gene mutations: a comparative clinical, histologic, and immunohistochemical study. *Arch. Dermatol.* 141: 343-351.

CHROMOSOMAL LOCATION

Genetic locus: HR (human) mapping to 8p21.3.

PRODUCT

HR siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see HR shRNA Plasmid (h): sc-77654-SH and HR shRNA (h) Lentiviral Particles: sc-77654-V as alternate gene silencing products.

For independent verification of HR (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-77654A, sc-77654B and sc-77654C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

HR siRNA (h) is recommended for the inhibition of HR expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

HR (A-11): sc-514686 is recommended as a control antibody for monitoring of HR gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor HR gene expression knockdown using RT-PCR Primer: HR (h)-PR: sc-77654-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.