

GalNAc-TL2 siRNA (h): sc-78461

BACKGROUND

The UDP-N-acetyl- α -D-galactosamine:polypeptide N-acetylgalactosaminyltransferase (GalNAc-T) family of enzymes are substrate-specific proteins that catalyze the transfer of GalNAc (N-acetylgalactosaminyl) to serine and threonine residues onto various proteins, thereby initiating mucin-type O-linked glycosylation in the Golgi apparatus. GalNAc-TL2 (polypeptide N-acetylgalactosaminyltransferase-like protein 2), also known as protein-UDP acetylgalactosaminyltransferase-like protein 2, is a 639 amino acid single-pass type II membrane protein that localizes to the Golgi apparatus and, like other GalNAc proteins, contains a stem region and a C-terminal ricin/lectin-like domain. GalNAc-TL2 utilizes calcium and manganese as cofactors. GalNAc-TL2 is widely expressed with highest expression in small intestine, placenta, spleen, cerebral cortex and ovary. Although it displays weaker catalytic activity than GalNAc-T2 and it prefers Muc1a as a substrate, GalNAc-TL2 can transfer up to seven GalNAc residues to the Muc5AC peptide.

REFERENCES

1. Porowska, H., Paszkiewicz-Gadek, A. and GindzieDski, A. 1999. Activity of partially purified UDP-N-acetyl- α -D-galactosamine: polypeptide N-acetylgalactosaminyltransferase with different peptide acceptors. *Acta Biochim. Pol.* 46: 365-370.
2. Kumar, S., Connor, J.R., Dodds, R.A., Halsey, W., Van Horn, M., Mao, J., Sathe, G., Mui, P., Agarwal, P., Badger, A.M., Lee, J.C., Gowen, M. and Lark, M.W. 2001. Identification and initial characterization of 5000 expressed sequenced tags (ESTs) each from adult human normal and osteoarthritic cartilage cDNA libraries. *Osteoarthr. Cartil.* 9: 641-653.
3. Schwientek, T., Bennett, E.P., Flores, C., Thacker, J., Hollmann, M., Reis, C.A., Behrens, J., Mandel, U., Keck, B., Schafer, M.A., Haselmann, K., Zubarev, R., Roepstorff, P., Burchell, J.M., Taylor-Papadimitriou, J., et al. 2002. Functional conservation of subfamilies of putative UDP-N-acetylgalactosamine:polypeptide N-acetylgalactosaminyltransferases in *Drosophila*, *Caenorhabditis elegans*, and mammals. One subfamily composed of I(2)35Aa is essential in *Drosophila*. *J. Biol. Chem.* 277: 22623-22638.
4. Arg_eso, P., Tisdale, A., Mandel, U., Letko, E., Foster, C.S. and Gipson, I.K. 2003. The cell-layer- and cell-type-specific distribution of GalNAc-transferases in the ocular surface epithelia is altered during keratinization. *Invest. Ophthalmol. Vis. Sci.* 44: 86-92.
5. Cheng, L., Tachibana, K., Iwasaki, H., Kameyama, A., Zhang, Y., Kubota, T., Hiruma, T., Tachibana, K., Kudo, T., Guo, J.M. and Narimatsu, H. 2004. Characterization of a novel human UDP-GalNAc transferase, pp-GalNAc-T15. *FEBS Lett.* 566: 17-24.
6. Herr, P., Korniyuchuk, G., Yamamoto, Y., Grubisic, K. and Oelgeschläger, M. 2008. Regulation of TGF- β signalling by N-acetylgalactosaminyltransferase-like 1. *Development* 135: 1813-1822.

CHROMOSOMAL LOCATION

Genetic locus: GALNTL2 (human) mapping to 3p25.1.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

GalNAc-TL2 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see GalNAc-TL2 shRNA Plasmid (h): sc-78461-SH and GalNAc-TL2 shRNA (h) Lentiviral Particles: sc-78461-V as alternate gene silencing products.

For independent verification of GalNAc-TL2 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-78461A, sc-78461B and sc-78461C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

GalNAc-TL2 siRNA (h) is recommended for the inhibition of GalNAc-TL2 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor GalNAc-TL2 gene expression knockdown using RT-PCR Primer: GalNAc-TL2 (h)-PR: sc-78461-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.