

GOS2 siRNA (h): sc-78689

BACKGROUND

GOS2 (putative lymphocyte G₀/G₁ switch protein 2) is a 103 amino acid novel target of peroxisome proliferator-activated receptors (PPARs) and regulator of latent HIV. GOS2 may be involved in adipocyte differentiation and its expression is essential for committing cells to enter the G₁ phase of the cell cycle. GOS2 contains a CpG-rich island and multiple sites for potential phosphorylation by casein kinase II and protein kinase C. The gene encoding GOS2 maps to human chromosome 1, which is the largest human chromosome. Chromosome 1 spans about 260 million base pairs and makes up 8% of the human genome. There are about 3,000 genes on chromosome 1 and, considering the great number of genes, there are also a large number of diseases associated with chromosome 1. Notably, the rare aging disease Hutchinson-Gilford progeria is associated with the LMNA gene which encodes lamin A. Stickler syndrome, Parkinson's, Gaucher disease and Usher syndrome are also associated with chromosome 1. Aberrations in chromosome 1 are found in a variety of cancers including head and neck cancer, malignant melanoma and multiple myeloma.

REFERENCES

1. Russell, L. and Forsdyke, D.R. 1991. A human putative lymphocyte G₀/G₁ switch gene containing a CpG-rich island encodes a small basic protein with the potential to be phosphorylated. *DNA Cell Biol.* 10: 581-591.
2. Cristillo, A.D., et al. 1997. Cyclosporin A inhibits early mRNA expression of G₀/G₁ switch gene 2 (GOS2) in cultured human blood mononuclear cells. *DNA Cell Biol.* 16: 1449-1458.
3. Eudy, J.D., et al. 1998. Mutation of a gene encoding a protein with extra-cellular matrix motifs in Usher syndrome type IIa. *Science* 280: 1753-1757.
4. Tayebi, N., et al. 2001. Gaucher disease and Parkinsonism: a phenotypic and genotypic characterization. *Mol. Genet. Metab.* 73: 313-321.
5. Oliveira, S.A., et al. 2005. Identification of risk and age-at-onset genes on chromosome 1p in Parkinson disease. *Am. J. Hum. Genet.* 77: 252-264.
6. Zandbergen, F., et al. 2005. The G₀/G₁ switch gene 2 is a novel PPAR target gene. *Biochem. J.* 392: 313-324.
7. Kitareewan, S., et al. 2008. GOS2 is an all-*trans*-retinoic acid target gene. *Int. J. Oncol.* 33: 397-404.
8. Balcáková, J., et al. 2009. Gain of chromosome arm 1q in patients in relapse and progression of multiple myeloma. *Cancer Genet. Cytogenet.* 192: 68-72.
9. Yokoi, T., et al. 2009. Analysis of the vitreous membrane in a case of type 1 Stickler syndrome. *Graefes Arch. Clin. Exp. Ophthalmol.* 247: 715-718.

CHROMOSOMAL LOCATION

Genetic locus: GOS2 (human) mapping to 1q32.2.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

GOS2 siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see GOS2 shRNA Plasmid (h): sc-78689-SH and GOS2 shRNA (h) Lentiviral Particles: sc-78689-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

GOS2 siRNA (h) is recommended for the inhibition of GOS2 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor GOS2 gene expression knockdown using RT-PCR Primer: GOS2 (h)-PR: sc-78689-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.