

WDTC1 siRNA (h): sc-88033

BACKGROUND

WDTC1 (WD and tetratricopeptide repeats 1), also known as ADP or DCAF9, is a 677 amino acid protein that contains seven WD repeats, two TPR repeats and exists as six alternatively spliced isoforms. WDTC1 is the human ortholog of *Drosophila* adp (adipose), an obesity gene involved in fat storage. WDTC1 is encoded by a gene located on human chromosome 1, which spans 260 million base pairs, contains over 3,000 genes and comprises nearly 8% of the human genome. Chromosome 1 houses a large number of disease-associated genes, including those involved in familial adenomatous polyposis, Stickler syndrome, Parkinson's disease, Gaucher disease, schizophrenia and Usher syndrome. Aberrations in chromosome 1 are found in a variety of cancers, including head and neck cancer, malignant melanoma and multiple myeloma.

REFERENCES

1. Eudy, J.D., Yao, S., Weston, M.D., Ma-Edmonds, M., Talmadge, C.B., Cheng, J.J., Kimberling, W.J. and Sumegi, J. 1998. Isolation of a gene encoding a novel member of the nuclear receptor superfamily from the critical region of Usher syndrome type IIa at 1q41. *Genomics* 50: 382-384.
2. Tayebi, N., Callahan, M., Madike, V., Stubblefield, B.K., Orvisky, E., Krasnewich, D., Fillano, J.J. and Sidransky, E. 2001. Gaucher disease and parkinsonism: a phenotypic and genotypic characterization. *Mol. Genet. Metab.* 73: 313-321.
3. Häder, T., Müller, S., Aguilera, M., Eulenberg, K.G., Steuernagel, A., Ciossek, T., Kühnlein, R.P., Lemaire, L., Fritsch, R., Dohrmann, C., Vetter, I.R., Jäckle, H., Doane, W.W. and Brönnner, G. 2003. Control of triglyceride storage by a WD40/TPR-domain protein. *EMBO Rep.* 4: 511-516.
4. Angers, S., Li, T., Yi, X., MacCoss, M.J., Moon, R.T. and Zheng, N. 2006. Molecular architecture and assembly of the DDB1-CUL4A ubiquitin ligase machinery. *Nature* 443: 590-593.
5. Lee, J. and Zhou, P. 2007. DCAFs, the missing link of the CUL4-DDB1 ubiquitin ligase. *Mol. Cell* 26: 775-780.
6. Balcárková, J., Urbánková, H., Scudla, V., Holzerová, M., Bacovský, J., Indrák, K. and Jarosová, M. 2009. Gain of chromosome arm 1q in patients in relapse and progression of multiple myeloma. *Cancer Genet. Cytogenet.* 192: 68-72.
7. Lai, C.Q., Parnell, L.D., Arnett, D.K., García-Bailo, B., Tsai, M.Y., Kabagambe, E.K., Straka, R.J., Province, M.A., An, P., Borecki, I.B., Tucker, K.L. and Ordovás, J.M. 2009. WDTC1, the ortholog of *Drosophila* adipose gene, associates with human obesity, modulated by MUFA intake. *Obesity* 17: 593-600.

CHROMOSOMAL LOCATION

Genetic locus: WDTC1 (human) mapping to 1p36.11.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

WDTC1 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see WDTC1 shRNA Plasmid (h): sc-88033-SH and WDTC1 shRNA (h) Lentiviral Particles: sc-88033-V as alternate gene silencing products.

For independent verification of WDTC1 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-88033A, sc-88033B and sc-88033C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

WDTC1 siRNA (h) is recommended for the inhibition of WDTC1 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor WDTC1 gene expression knockdown using RT-PCR Primer: WDTC1 (h)-PR: sc-88033-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.