

SCaMC-1 siRNA (h): sc-88209

BACKGROUND

SCaMC-1 (short calcium-binding mitochondrial carrier 1), also known as APC1 or SLC25A24 (solute carrier family 25, member 24), is a 477 amino acid multi-pass membrane protein belonging to the SLC25 family of mitochondrial carriers, which are responsible for transporting metabolites across the inner mitochondrial membrane. Expressed in a wide variety of tissues and localizing to the mitochondrial inner membrane, SCaMC-1 contains three Solcar repeats and four EF-hand domains and functions as a calcium-dependent mitochondrial solute carrier. SCaMC-1 may act as an ATP-Mg/Pi co-transporter, facilitating the transport of Mg-ATP in exchange for phosphate. Existing as two isoforms, SCaMC-1 is encoded by a gene located on human chromosome 1, which is the largest human chromosome spanning about 260 million base pairs and making up 8% of the human genome.

REFERENCES

1. Mashima, H., et al. 2003. A novel mitochondrial Ca^{2+} -dependent solute carrier in the liver identified by mRNA differential display. *J. Biol. Chem.* 278: 9520-9527.
2. del Arco, A., et al. 2004. Identification of a novel human subfamily of mitochondrial carriers with calcium-binding domains. *J. Biol. Chem.* 279: 24701-24713.
3. Fiermonte, G., et al. 2004. Identification of the mitochondrial ATP-Mg/Pi transporter. Bacterial expression, reconstitution, functional characterization, and tissue distribution. *J. Biol. Chem.* 279: 30722-30730.
4. Online Mendelian Inheritance in Man, OMIM™. 2004. Johns Hopkins University, Baltimore, MD. MIM Number: 608744. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Del Arco, A. 2005. Novel variants of human SCaMC-3, an isoform of the ATP-Mg/Pi mitochondrial carrier, generated by alternative splicing from 3'-flanking transposable elements. *Biochem. J.* 389: 647-655.
6. Bassi, M.T., et al. 2005. Cellular expression and alternative splicing of SLC25A23, a member of the mitochondrial Ca^{2+} -dependent solute carrier gene family. *Gene* 345: 173-182.
7. Gregory, S.G., et al. 2006. The DNA sequence and biological annotation of human chromosome 1. *Nature* 441: 315-321.

CHROMOSOMAL LOCATION

Genetic locus: SLC25A24 (human) mapping to 1p13.3.

PRODUCT

SCaMC-1 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μM solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see SCaMC-1 shRNA Plasmid (h): sc-88209-SH and SCaMC-1 shRNA (h) Lentiviral Particles: sc-88209-V as alternate gene silencing products.

For independent verification of SCaMC-1 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-88209A, sc-88209B and sc-88209C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20°C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20°C , avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μl of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μl of RNase-free water makes a 10 μM solution in a 10 μM Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

SCaMC-1 siRNA (h) is recommended for the inhibition of SCaMC-1 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μM in 66 μl . Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor SCaMC-1 gene expression knockdown using RT-PCR Primer: SCaMC-1 (h)-PR: sc-88209-PR (20 μl). Annealing temperature for the primers should be $55-60^{\circ}\text{C}$ and the extension temperature should be $68-72^{\circ}\text{C}$.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.