

# DHDDS siRNA (h): sc-88357

## BACKGROUND

DHDDS (dehydrodolichyl diphosphate synthase), also known as CIT, CPT, HDS or Dedol-PP synthase, is a 333 amino acid protein that catalyzes *cis*-prenyl chain elongation to form the polyprenyl backbone of dolichol. A member of the UPP synthase family, DHDDS is highly expressed in kidney and testis, with lower levels found in spleen, thymus and heart. DHDDS exists as three alternatively spliced isoforms and is encoded by a gene that maps to human chromosome 1p36.11. Human chromosome 1 spans 260 million base pairs, contains over 3,000 genes, comprises nearly 8% of the human genome and houses a large number of disease-associated genes, including those that are involved in familial adenomatous polyposis, Stickler syndrome, Parkinson's disease, Gaucher disease, schizophrenia and Usher syndrome.

## REFERENCES

1. Eudy, J.D., Weston, M.D., Yao, S., Hoover, D.M., Rehm, H.L., Ma-Edmonds, M., Yan, D., Ahmad, I., Cheng, J.J., Ayuso, C., Cremers, C., Davenport, S., Moller, C., Talmadge, C.B., Beisel, K.W., Tamayo, M., et al. 1998. Mutation of a gene encoding a protein with extracellular matrix motifs in Usher syndrome type IIa. *Science* 280: 1753-1757.
2. Tayebi, N., Callahan, M., Madike, V., Stubblefield, B.K., Orvisky, E., Krasnewich, D., Fillano, J.J. and Sidransky, E. 2001. Gaucher disease and parkinsonism: a phenotypic and genotypic characterization. *Mol. Genet. Metab.* 73: 313-321.
3. Endo, S., Zhang, Y.W., Takahashi, S. and Koyama, T. 2003. Identification of human dehydrodolichyl diphosphate synthase gene. *Biochim. Biophys. Acta* 1625: 291-295.
4. Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2003. Johns Hopkins University, Baltimore, MD. MIM Number: 608172. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Pilaslova, M., Russell, A.M., Wanner, A., Wolf, A., Dobbie, Z., Müller, H.J. and Heinemann, K. 2004. Exclusion of an extracolonic disease modifier locus on chromosome 1p33-36 in a large Swiss familial adenomatous polyposis kindred. *Eur. J. Hum. Genet.* 12: 365-371.
6. Jones, J., Viswanathan, K., Krag, S.S. and Betenbaugh, M.J. 2005. Polyprenyl lipid synthesis in mammalian cells expressing human *cis*-prenyl transferase. *Biochem. Biophys. Res. Commun.* 331: 379-383.
7. Yurov, Y.B., Iourov, I.Y., Vorsanova, S.G., Demidova, I.A., Kravetz, V.S., Beresheva, A.K., Kolotii, A.D., Monakhov, V.V., Uranova, N.A., Vostrikov, V.M., Soloviev, I.V. and Liehr, T. 2008. The schizophrenia brain exhibits low-level aneuploidy involving chromosome 1. *Schizophr. Res.* 98: 139-147.
8. Yokoi, T., Koide, R., Matsuoka, K., Nakagawa, A. and Azuma, N. 2009. Analysis of the vitreous membrane in a case of type 1 Stickler syndrome. *Graefes Arch. Clin. Exp. Ophthalmol.* 247: 715-718.
9. Li, J., Liu, F., Wang, H., Liu, X., Liu, J., Li, N., Wan, F., Wang, W., Zhang, C., Jin, S., Liu, J., Zhu, P. and Liu, Y. 2010. Systematic mapping and functional analysis of a family of human epididymal secretory sperm-located proteins. *Mol. Cell. Proteomics* 9: 2517-2528.

## CHROMOSOMAL LOCATION

Genetic locus: DHDDS (human) mapping to 1p36.11.

## PRODUCT

DHDDS siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see DHDDS shRNA Plasmid (h): sc-88357-SH and DHDDS shRNA (h) Lentiviral Particles: sc-88357-V as alternate gene silencing products.

For independent verification of DHDDS (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-88357A, sc-88357B and sc-88357C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

DHDDS siRNA (h) is recommended for the inhibition of DHDDS expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor DHDDS gene expression knockdown using RT-PCR Primer: DHDDS (h)-PR: sc-88357-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.