

# 17 $\beta$ -HSD7 siRNA (h): sc-88433

## BACKGROUND

17 $\beta$ -HSD7 (17 $\beta$  hydroxysteroid dehydrogenase type 7), also designated 3-ketosteroid reductase, belongs to the 17 $\beta$ -HSD family of proteins, which regulate the availability of steroids within various tissues throughout the body. 17 $\beta$ -HSD7 is a 341 amino acid protein that converts estrone to estradiol and is also involved in cholesterol biosynthesis. 17 $\beta$ -HSD7 is highly expressed in adrenal gland, liver, lung and thymus. It is also expressed in the corpus luteum, where it is thought to play a role in fetal development. Single nucleotide polymorphisms in the gene encoding 17 $\beta$ -HSD7 have been shown to affect its level of transcription in LNCaP and DU145 cells, which may modulate an adverse reaction induced by estramustine phosphate sodium.

## REFERENCES

1. Nokelainen, P., Peltoketo, H., Vihko, R. and Vihko, P. 1998. Expression cloning of a novel estrogenic mouse 17 $\beta$ -hydroxysteroid dehydrogenase/17-ketosteroid reductase (m17HSD7), previously described as a Prolactin receptor-associated protein (PRAP) in rat. *Mol. Endocrinol.* 12: 1048-1059.
2. Marijanovic, Z., Laubner, D., Moller, G., Gege, C., Husen, B., Adamski, J. and Breitling, R. 2003. Closing the gap: identification of human 3-ketosteroid reductase, the last unknown enzyme of mammalian cholesterol biosynthesis. *Mol. Endocrinol.* 17: 1715-1725.
3. Risk, M., Shehu, A., Mao, J., Stocco, C.O., Goldsmith, L.T., Bowen-Shauver, J.M. and Gibori, G. 2005. Cloning and characterization of a 5' regulatory region of the Prolactin receptor-associated protein/17 $\beta$  hydroxysteroid dehydrogenase 7 gene. *Endocrinology* 146: 2807-2816.
4. Seth, G., Mclvor, R.S. and Hu, W.S. 2006. 17 $\beta$ -hydroxysteroid dehydrogenase type 7 (Hsd17b7) reverts cholesterol auxotrophy in NS0 cells. *J. Biotechnol.* 121: 241-252.
5. Ohnesorg, T., Keller, B., Hrab $\acute{e}$  de Angelis, M. and Adamski, J. 2006. Transcriptional regulation of human and murine 17 $\beta$ -hydroxysteroid dehydrogenase type-7 confers its participation in cholesterol biosynthesis. *J. Mol. Endocrinol.* 37: 185-197.

## CHROMOSOMAL LOCATION

Genetic locus: HSD17B7 (human) mapping to 1q23.3.

## PRODUCT

17 $\beta$ -HSD7 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see 17 $\beta$ -HSD7 shRNA Plasmid (h): sc-88433-SH and 17 $\beta$ -HSD7 shRNA (h) Lentiviral Particles: sc-88433-V as alternate gene silencing products.

For independent verification of 17 $\beta$ -HSD7 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-88433A, sc-88433B and sc-88433C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20 $^{\circ}$  C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20 $^{\circ}$  C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

17 $\beta$ -HSD7 siRNA (h) is recommended for the inhibition of 17 $\beta$ -HSD7 expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

17 $\beta$ -HSD7 (F-4): sc-393936 is recommended as a control antibody for monitoring of 17 $\beta$ -HSD7 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker<sup>TM</sup> Molecular Weight Standards: sc-2035, UltraCruz<sup>®</sup> Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz<sup>®</sup> Mounting Medium: sc-24941 or UltraCruz<sup>®</sup> Hard-set Mounting Medium: sc-359850.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor 17 $\beta$ -HSD7 gene expression knockdown using RT-PCR Primer: 17 $\beta$ -HSD7 (h)-PR: sc-88433-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60 $^{\circ}$  C and the extension temperature should be 68-72 $^{\circ}$  C.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.