

# CYP4Z1 siRNA (h): sc-88518

## BACKGROUND

Cytochrome P450 proteins are heme-thiolate monooxygenases that mediate NADPH-dependent electron transport and function to oxidize a variety of structurally unrelated compounds, including steroids, fatty acids and xenobiotics. Specifically, cytochrome P450s are responsible for metabolizing arachidonic acid to hydroxyeicosatetraenoic acid (a regulator of blood pressure) and epoxyeicosatrienoic acid (a molecule involved in signaling events). CYP4Z1 (cytochrome P450, family 4, subfamily Z, polypeptide 1), also known as CYP4A20, is a 505 amino acid single-pass type II membrane protein that localizes to both the microsome and the endoplasmic reticulum and belongs to the cytochrome P450 family. Using heme groups as cofactors, CYP4Z1 functions to catalyze the conversion of a reduced flavoprotein to an oxidized flavoprotein and is involved in drug metabolism, as well as in the synthesis of cholesterol, steroids and other lipids. CYP4Z1 is expressed at high levels in breast cancer tissue, suggesting a role in tumorigenesis.

## REFERENCES

1. Simpson, A.E. 1997. The cytochrome P450 4 (CYP4) family. *Gen. Pharmacol.* 28: 351-359.
2. Rieger, M.A., Ebner, R., Bell, D.R., Kiessling, A., Rohayem, J., Schmitz, M., Temme, A., Rieber, E.P. and Weigle, B. 2004. Identification of a novel mammary-restricted cytochrome P450, CYP4Z1, with overexpression in breast carcinoma. *Cancer Res.* 64: 2357-2364.
3. Nelson, D.R., Zeldin, D.C., Hoffman, S.M., Maltais, L.J., Wain, H.M. and Nebert, D.W. 2004. Comparison of cytochrome P450 (CYP) genes from the mouse and human genomes, including nomenclature recommendations for genes, pseudogenes and alternative-splice variants. *Pharmacogenetics* 14: 1-18.
4. Savas, U., Hsu, M.H., Griffin, K.J., Bell, D.R. and Johnson, E.F. 2005. Conditional regulation of the human CYP4X1 and CYP4Z1 genes. *Arch. Biochem. Biophys.* 436: 377-385.
5. Thomas, R.D., Green, M.R., Wilson, C., Weckle, A.L., Duanmu, Z., Kocarek, T.A. and Runge-Morris, M. 2006. Cytochrome P450 expression and metabolic activation of cooked food mutagen 2-amino-1-methyl-6-phenylimidazo[4,5-b]pyridine (PhIP) in MCF10A breast epithelial cells. *Chem. Biol. Interact.* 160: 204-216.
6. Zöllner, A., Dragan, C.A., Pistorius, D., Müller, R., Bode, H.B., Peters, F.T., Maurer, H.H. and Bureik, M. 2009. Human CYP4Z1 catalyzes the in-chain hydroxylation of lauric acid and myristic acid. *Biol. Chem.* 390: 313-317.
7. Bankovic, J., Stojic, J., Jovanovic, D., Andjelkovic, T., Milinkovic, V., Ruzdijic, S. and Tanic, N. 2010. Identification of genes associated with non-small-cell lung cancer promotion and progression. *Lung Cancer* 67: 151-159.

## CHROMOSOMAL LOCATION

Genetic locus: CYP4Z1 (human) mapping to 1p33.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

## PRODUCT

CYP4Z1 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see CYP4Z1 shRNA Plasmid (h): sc-88518-SH and CYP4Z1 shRNA (h) Lentiviral Particles: sc-88518-V as alternate gene silencing products.

For independent verification of CYP4Z1 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-88518A, sc-88518B and sc-88518C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

CYP4Z1 siRNA (h) is recommended for the inhibition of CYP4Z1 expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor CYP4Z1 gene expression knockdown using RT-PCR Primer: CYP4Z1 (h)-PR: sc-88518-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.