

SBZF3 siRNA (h): sc-88521

BACKGROUND

Zinc-finger proteins contain DNA-binding domains and have a wide variety of functions, most of which encompass some form of transcriptional activation or repression. The majority of zinc-finger proteins contain a Krüppel-type DNA binding domain and a KRAB domain, which is thought to interact with KAP1, thereby recruiting histone modifying proteins. SBZF3, also known as ZNF695 (zinc finger protein 695), is a 515 amino acid nuclear protein that contains 13 C₂H₂-type zinc fingers and one KRAB domain. Like other zinc finger proteins, it is thought to play a role in transcriptional regulation. The gene encoding SBZF3 maps to human chromosome 1q44, the largest human chromosome spanning about 260 million base pairs and making up 8% of the human genome. There are four isoforms of SBZF3 that are produced as a result of alternative splicing events.

REFERENCES

1. Matthews, J.M. and Sunde, M. 2002. Zinc fingers—folds for many occasions. *IUBMB Life* 54: 351-355.
2. Hurley, J.H., Lee, S. and Prag, G. 2006. Ubiquitin-binding domains. *Biochem. J.* 399: 361-372.
3. Gamsjaeger, R., Liew, C.K., Loughlin, F.E., Crossley, M. and Mackay, J.P. 2007. Sticky fingers: zinc-fingers as protein-recognition motifs. *Trends Biochem. Sci.* 32: 63-70.
4. Itokawa, Y., Yanagawa, T., Yamakawa, H., Watanabe, N., Koga, H. and Nagase, T. 2009. KAP1-independent transcriptional repression of SCAN-KRAB-containing zinc finger proteins. *Biochem. Biophys. Res. Commun.* 388: 689-694.
5. Ding, G., Lorenz, P., Kreutzer, M., Li, Y. and Thiesen, H.J. 2009. SysZNF: the C₂H₂ zinc finger gene database. *Nucleic Acids Res.* 37: D267-D273.
6. Lorenz, P., Dietmann, S., Wilhelm, T., Koczan, D., Autran, S., Gad, S., Wen, G., Ding, G., Li, Y., Rousseau-Merck, M.F. and Thiesen, H.J. 2010. The ancient mammalian KRAB zinc finger gene cluster on human chromosome 8q24.3 illustrates principles of C₂H₂ zinc finger evolution associated with unique expression profiles in human tissues. *BMC Genomics* 11: 206.
7. Fietze, S., Lan, X., Jin, V.X. and Farnham, P.J. 2010. Genomic targets of the KRAB and SCAN domain-containing zinc finger protein 263. *J. Biol. Chem.* 285: 1393-1403.
8. Nowick, K., Hamilton, A.T., Zhang, H. and Stubbs, L. 2010. Rapid sequence and expression divergence suggests selection for novel function in primate-specific KRAB-ZNF genes. *Mol. Biol. Evol.* 27: 2606-2617.

CHROMOSOMAL LOCATION

Genetic locus: ZNF695 (human) mapping to 1q44.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

SBZF3 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see SBZF3 shRNA Plasmid (h): sc-88521-SH and SBZF3 shRNA (h) Lentiviral Particles: sc-88521-V as alternate gene silencing products.

For independent verification of SBZF3 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-88521A, sc-88521B and sc-88521C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

SBZF3 siRNA (h) is recommended for the inhibition of SBZF3 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor SBZF3 gene expression knockdown using RT-PCR Primer: SBZF3 (h)-PR: sc-88521-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.