

p52 S6 kinase siRNA (h): sc-88609

BACKGROUND

The phosphorylation and dephosphorylation of proteins on serine and threonine residues is an essential means of regulating a broad range of cellular functions in eukaryotes, including cell division, homeostasis and apoptosis. A group of proteins that are intimately involved in this process are the serine/threonine (Ser/Thr) protein kinases. p52 S6 kinase, also known as RPS6KC1 (ribosomal protein S6 kinase, 52 kDa, polypeptide 1) or RPK118, is a 1,066 amino acid member of the Ser/Thr kinase family that localizes to both the cytoplasm and the nucleus and contains one MIT domain, one PX domain and two protein kinase domains. Expressed at high levels in brain, placenta, heart, testis, kidney, liver and skeletal muscle, p52 S6 kinase catalyzes the ATP-dependent phosphorylation of target proteins and is thought to be involved in transmitting sphingosine-1 phosphate (SPP)-mediated signaling into the cell.

REFERENCES

- Loijens, J.C., et al. 1996. The phosphatidylinositol 4-phosphate 5-kinase family. *Adv. Enzyme Regul.* 36: 115-140.
- Zhang, H., et al. 1999. Cloning, characterization, and chromosome mapping of RPS6KC1, a novel putative member of the ribosome protein S6 kinase family, to chromosome 12q12-q13.1. *Genomics* 61: 314-318.
- Hayashi, S., et al. 2002. Identification and characterization of RPK118, a novel sphingosine kinase-1-binding protein. *J. Biol. Chem.* 277: 33319-33324.
- Ellson, C.D., et al. 2002. The PX domain: a new phosphoinositide-binding module. *J. Cell Sci.* 115: 1099-1105.
- Liu, H., et al. 2002. Sphingosine kinases: a novel family of lipid kinases. *Prog. Nucleic Acid Res. Mol. Biol.* 71: 493-511.
- Ishida, S., et al. 2004. Differential modulation of PI 3-kinase/Akt pathway during all-trans retinoic acid- and Am80-induced HL-60 cell differentiation revealed by DNA microarray analysis. *Biochem. Pharmacol.* 68: 2177-2186.
- Liu, L., et al. 2005. RPK118, a PX domain-containing protein, interacts with peroxiredoxin-3 through pseudo-kinase domains. *Mol. Cells* 19: 39-45.

CHROMOSOMAL LOCATION

Genetic locus: RPS6KC1 (human) mapping to 1q32.3.

PRODUCT

p52 S6 kinase siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see p52 S6 kinase shRNA Plasmid (h): sc-88609-SH and p52 S6 kinase shRNA (h) Lentiviral Particles: sc-88609-V as alternate gene silencing products.

For independent verification of p52 S6 kinase (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-88609A, sc-88609B and sc-88609C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

p52 S6 kinase siRNA (h) is recommended for the inhibition of p52 S6 kinase expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor p52 S6 kinase gene expression knockdown using RT-PCR Primer: p52 S6 kinase (h)-PR: sc-88609-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.